

**Q1: How long is the power cord? Does it plug into the wall socket?**

The power cord of this product is roughly 5ft long and it can plug into the wall socket.

**Q2: Is it desktop only?**

Not at all. There is a hanging hole on the back of the product, which can be optionally placed on the table or hung on the wall.

**Q3: Will the light from this light tan my face?**

No, there is no UV hazard associated with the use of our products.

**Q4: Does it need to be plugged in all the time?**

Yes, the light therapy lamp must be plugged in to work. This item comes with a power adapter, and the power supply is designed to accept anywhere from 100V to 240V.

**Q5: Does it have a timer so I can set it to turn on and then off daily at specific times?**

No, it does not have this feature, but the device has a timer for setting use-time and will shut down when time is over.

**Q6:What is sunset mode?**

Sunset mode is holding B and C simultaneously. After choosing soft sunset mode, the lights dims over 20 minutes period until it slowly fades out, perfectly mirroring an evening sunset.

**Q7: Is it made of wood?**

No. It's made of acrylic, but the lamp has exquisite faux wood grain.

We hope these tips were useful. If you have any inquiries about our products or services feel free to contact us via email, call our customer service at 855.822.6999 M-F 9am-5pm CT or visit us at [healthcare-manager.com](https://www.healthcare-manager.com)