

WUJIANG PAILIXIN METAL PRODUCT CO; LTD

Sterilization Certificate

General Information

Sterilization No.	190201-1		
Customer name	PEAK		
Date of reception	2019.02.01		
Temperature of Environment	50°C		
Product name	Product type	Product Lot#	Cases
STAINLESS STEEL PIN WITH TUBE	1RL 0.18MM、1RL 0.25MM、 1RL 0.35MM、1RL 0.40MM etc	0201	32CTNS

Sterilization process data

Preconditioning time	1.5-2.5hrs
Date of sterilization	2019.02.01
Station number	7
Batch number	190201
Aeration time	7 Days

Results of Physical process control

The physical Acceptance criteria of the preconditioning the cycl and the aeration were in compliance with the validated specifications of cycle no.:GBS

Microbiological process data

The microbiological process control with biological indicators has been performed according to the BI sterility test sop(032005)

No of biological indicators indicators 19 per chamber load.	Beijin Four Ring Medical, Lot No.20170203
Date of incubation	2019.02.01
Incubation conditions	7 Days 37°C

Result of microbiological process control

[] No growth from Bacillus subtilis has been detected

Remarks

--	--

QA Approval

Date of release	Signature of QA
2019.02.01	