

ELECTRIC WHEELCHAIRS



ULTRA-LITE 2 ELECTRIC WHEELCHAIR

- LTA compliant personal mobility aid (PMA).
- Each wheelchair comes with 1 quick-release lithium battery. Tools not required.
- Lightweight, only 16 kg without battery and footrest. Fully loaded weight only 18 kg.
- 24V, 6 AH lithium battery that travel up to 10 km.
- Top speed of 6 km/h
- Flip-up armrests and detachable footrests.
- Front wheel suspension.

SPECIFICATIONS

Weight Capacity	: 100 kg
Frame Weight	: 16.1 kg (w/o battery & footrests)
Battery Weight	: 1.3 kg x 1 pieces
Footrests Weight	: 0.5 kg
Total Weight (Fully Loaded)	: 17.9 kg
Overall Width	: 55 cm
Overall Length	: 91 cm
Overall Height	: 89 cm
Minimum Folded Dimensions	: 71.5 cm x 37 cm x 73.5 cm
Seat Size (W x D x H)	: 44 cm x 39 cm x 40 cm
Armrest Height	: 23 cm
Seat-to-Floor Height	: 47 cm
Travelling Range	: 10 km per battery
Max Speed	: 6 km/h
Turning Radius	: 26"
Batteries	: 1 x 24V, 6 AH Lithium
Charger	: 2A c/w Safety Mark
Motors	: 24V 150W x 2
Max Gradient	: 6 degrees
Wheels	: Front 6" Solid, Rear 12" Solid