

# MOBILITY SCOOTERS

## BUDGET-LITE DELUXE LARGE 4-WHEELED MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Captain's seat that can swivel 360 degrees
- Backrest angle, headrest height and tiller angle can be adjusted anytime without tools
- Flip-up Armrests
- Full suspension system with large tyres and high ground clearance
- Full lighting system with front, rear, signal and hazard lights
- Rear view mirrors included



### SPECIFICATIONS

Overall Dimensions (mm)	: 1200 (L) x 600 (W) x 1013 (H)
Gross Weight with Battery	: 82 kg
Turning Radius	: 1200 mm
Max Speed	: 6 km/h
Max Climbing Angle	: 12 degrees
Max Range	: 30 km
Max Load	: 136 kg
Motor	: 400W / 24V
Battery	: 2 x 12V, 38AH
Charger	: 6A c/w Safety Mark
Front Wheels	: 8.5" Solid
Rear Wheels	: 10" Solid
Ground Clearance	: 80 mm
Seat Width	: 17"
Seat Depth	: 18"
Backrest Height	: 27" - 32" Adjustable ((including headrest)
Seat-to-Deck Height	: 14.5" - 17.5" Adjustable (tools required)
Seat-to-Floor Height	: 20" - 23" Adjustable (tools required)