

# MOBILITY SCOOTERS



## SW1000S PMA MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Adjustable swivel seat height and armrest width
- Quick detachable battery box
- Off-board charging
- Maintenance-free, puncture-proof solid tyres

## SPECIFICATIONS

Size (Unfolded)	: 112 cm x 54 cm x 96 cm
Base Weight	: 24.3 kg
Seat Weight	: 6.5 kg
Battery Weight	: 8.9 kg (12AH) / 13.9 kg (20AH)
Basket Weight	: 0.5 kg
Total Weight	: 40.2 kg (12AH) / 45.2 kg (20AH)
Max Capacity	: 120 kg
Max Speed	: 8 km/h
Max Range	: 12 km (12AH) / 20 km (20AH)
Max Gradient	: 10 degrees
Turning Radius	: 150 cm
Ground Clearance	: 4.5 cm
Seat Size	: 43 cm (width) x 38 cm (depth)
Backrest Height	: 33 cm
Battery	: 2 x 12V, 12AH / 2 x 12V, 20AH SLA
Motor	: 300w
Charger	: 24V 2A
Braking System	: Electromagnetic
Front Wheels	: 215 mm x 50 mm Solid
Rear Wheels	: 215 mm x 90 mm Solid