

MANUAL WHEELCHAIR



4-WHEEL ALUMINIUM ROLLATOR

- Handle height can be adjusted in 1" increments.
- Comes with new seamless padded seat with zippered pouch under seat.
- Removable, hinged padded backrest can be folded up or down.
- Comes with standard carry pouch under the seat for added privacy and security of personal items.
- Easy to use deluxe loop locks.
- Serrated brakes.
- Ergonomic handles are easy to grip and are adjustable in height.

SPECIFICATIONS

Overall Width	: 61 cm
Seat Width	: 45 cm
Rear Wheel	: 15 cm
Front Wheel	: 15 cm
Seat Height	: 55 cm
Overall Length	: 61 cm
Backrest Height	: 25 cm
Overall Height	: 79 - 90 cm adjustable
Folded Dimensions	: 63 x 24 x 85 cm
Loading Capacity	: 100 kg
Net Weight	: 6.9 kg