

MANUAL WHEELCHAIR

ALUMINIUM PUSHCHAIR

- Lightweight
- Easily foldable
- Comes with a fold-down backrest that makes it compact enough to fit most car boots
- Assisted brakes that allow caregivers to safely push the user around
- Rear wheels are made of high quality solid Polyurethane (PU) tyres, which are more wear resistant and provide better shock absorption than conventional solid tyres



SPECIFICATIONS

Overall Width	: 59 cm
Seat Size (Width x Depth)	: 46 cm x 38 cm
Rear Wheel	: 31 cm, solid PU
Front Wheel	: 15 cm, solid
Seat Height	: 49 cm
Overall Height	: 90 cm
Folded Dimensions	: 71 x 25 x 71 cm
Loading Capacity	: 100 kg
Net Weight	: 9.6 kg