LC0401S HepG2 Cytosol

Cytosolic lysate boiled in SDS-PAGE sample buffer prepared from HepG2 cells that were cultivated to 80% confluency and deprived of serum for 18-20 hours prior to harvesting



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

prior 20 2
,

Applications

Lysate Use Description: For testing antibodies by immunoblotting.

This product is for in vitro research use only and is not intended for use in humans or animals.