



Check out all our products at:
www.GoldenWestDiagnostics.com

- ◆ *Mass Spect Gold Human Matrices*
- ◆ *Processed Human Serums*
- ◆ *Protein Fractions*
- ◆ *OEM/Contract Manufacturing*



**ANIMAL PROTEIN
FRACTIONS**

Golden West Diagnostics, LLC. (GWD) – est. 2015. As a subsidiary of Golden West Biologicals, Inc (GWB), GWD is rooted in over 30 years of industry experience. GWD specializes in genuine human matrices used in the preparation for standards and controls for mass spectrometry and immunoassay applications. At GWD, we offer solutions to meet customers unique requirements. Our products are of the highest quality.

951-414-1870 | Orders@GoldenWestDiagnostics.com

ANIMAL PROTEIN FRACTIONS

Golden West Diagnostics (GWD) consistently supplies high-quality animal protein fractions. GWD offers ready to buy and custom products. Call us today with your specifications.

ANIMAL GAMMA GLOBULIN PRODUCTS

AGG2030 GOAT GAMMA GLOBULIN

TEST	SPECIFICATIONS	METHODOLOGY
Purity*	≥95%	HPLC
Total Protein*	≥4.5 g/dL	Biuret
pH*	6.0 - 8.0	pH Meter
Sodium*	50-200 mmol/L	ISE
Potassium*	<1.0 mmol/L	ISE
Chloride*	30-200 mmol/L	ISE
Moisture*	≤5%	Protein Calculation
Microbial*	<100 CFU/mL	USP/BAM
Calcium*	≤2 mg/dL	Arsenazo III

*5% Solution

AGG2090 RABBIT GAMMA GLOBULIN

TEST	SPECIFICATIONS	METHODOLOGY
Purity*	≥95%	HPLC
Total Protein*	≥4.5 g/dL	Biuret
pH*	6.0 - 8.0	pH Meter
Sodium*	50-200 mmol/L	ISE
Potassium*	<1.0 mmol/L	ISE
Chloride*	30-200 mmol/L	ISE
Moisture*	≤5%	Protein Calculation
Microbial*	<100 CFU/mL	USP/BAM
Calcium*	≤2 mg/dL	Arsenazo III

*5% Solution

AGG2110 SHEEP GAMMA GLOBULIN

TEST	SPECIFICATIONS	METHODOLOGY
Purity*	≥95%	HPLC
Total Protein*	≥4.5 g/dL	Biuret
pH*	6.0 - 8.0	pH Meter
Sodium*	50-200 mmol/L	ISE
Potassium*	<1.0 mmol/L	ISE
Chloride*	30-200 mmol/L	ISE
Moisture*	≤5%	Protein Calculation
Microbial*	<100 CFU/mL	USP/BAM
Calcium*	≤2 mg/dL	Arsenazo III

*5% Solution

THROMBIN

BF3040 BOVINE THROMBIN

Total Activity Per Container >1,000,000 US Unit

