

Disposable Serological Pipettes



PP Dispo Pipette



Dispo Pipette



PP Dispo Pipette

PP Dispo Pipette made of polypropylene is characterized by high resistance to chemicals and solvents and a low cost. A wider range of applications!

Best chemical resistance

Can be used for various purposes.

Low Cost

PP Dispo Pipette is less expensive because of polypropylene, compared with dispo pipette made of polystyrene or glass.

***After autoclaving this product, we CANNOT guarantee scale accuracy.**



High safety

PP Dispo Pipette does not break easily.


Short length

Easy to handle on a clean bench ensuring efficient operation.

- 5 mL Length 197 mm (With adapter length 220 mm)
- 10 mL Length 216.9 mm (With adapter length 239 mm)



***WATSON PP Dispo Pipette is also applicable in other brands' dispensers.**



Each adaptor is made of PP and contains PE filter.



Material: Polypropylene Body only

Cat. No.	Item	Unit
1653-005	5 mL PP Dispo Pipette with Dedicated adapter (1 piece)	100 pieces × 5 bags
1653-010	10 mL PP Dispo Pipette with Dedicated adapter (1 piece)	100 pieces × 3 bags
1653-AD5	Dedicated adapter for 5 mL PP Dispo Pipette	1 piece
1653-AD10	Dedicated adapter for 10 mL PP Dispo Pipette	1 piece

PP Dispo Pipette / Dispo Pipette

Dispo Pipette

High precision scale realized by our injection molding technology.

High precision scales



Precise & easily-recognizable scale realized by our advanced injection molding technology.

2 kinds of scales



Easy to see aspirated or dispensed volume with 2 different kinds of scales employed.

5 mL
Length : 224 mm

10 mL
Length : 244 mm



Each adaptor is made of PP and contains PE filter.

Short length

Easy-to-handle on a clean bench with limited handling space.
Length: 224 mm for 5 mL and 244 mm for 10 mL.

Clean

Sterilized, RNase/DNase-free
Endotoxin level: less than 0.05 EU/mL
Source: Internal inspection.



Material: Polystyrene (body) Polypropylene (adapter) Polyethylene (filter)

Cat. No.	Item	Unit	Case
1650-005S	5 mL Disposable Serological Pipette Sterilized	150 pieces (individual packaging)	5 units
1650-010S	10 mL Disposable Serological Pipette Sterilized	100 pieces (individual packaging)	5 units



Space-saving and easy-to-stack package design.