



# Racks/Stand/Cases



for Tubes




for Microscope Slide

Tube Stands



Material: Polypropylene(100% recycled)




The bottom of each well is shaped to stabilize the tube.

ECO Rack



- Ecology and economy.
- Made of 100% recycled polypropylene material.
- 96-well rack for 1.5 / 2.0 mL tubes, one size larger than 80 well rack.
- Writing space available on the side of the rack.



96-well rack with each hole obliquely-lined design.



18 mm pitch well interval. Multi channel pipette that installs the tip by skipping by one channel can be used.

Cat. No.	Item	Color	Dimensions	Case
1521-902B	Microtube Stand for 96 tubes "ECO RACK"	Blue	240.2 × 76.7 × 26.5 mm	50 pieces
1521-902BL	Microtube Stand for 96 tubes "ECO RACK"	Black	240.2 × 76.7 × 26.5 mm	50 pieces
1521-902W	Microtube Stand for 96 tubes "ECO RACK"	Natural	240.2 × 76.7 × 26.5 mm	50 pieces
1521-902Y	Microtube Stand for 96 tubes "ECO RACK"	Yellow	240.2 × 76.7 × 26.5 mm	50 pieces

Clear Tube Stand



Material: Polycarbonate



- Enable to set 20 pcs of Microtubes (only 1.5 mL applicable).
- The design is 2 rows with stair tyle.
- The bottom of tube is clearly visible while setting up.
- Anti-slip rubber equipped in two corners.



\*:see below

Cat. No.	Item	Dimensions	Unit	Case
1511-205	Microtube Stand for 20 tubes, Clear	210 × 65 × 35 mm	10 pieces	5 units


Note:(\*)Repeated autoclaving may cause rack damaged.



Piccolo Racks




1522-2451



1522-2453

- Mini tube rack for 24 tubes.
- Autoclavable rack made of polycarbonate.
- Possible to store in the liquid nitrogen (vapor phase) and the deep freezing condition.



The bottom of the tube is clearly visible while setting up.  
It is possible to set the crashed ice into the empty space of the rack, then keep the sample cool condition.



The notch(Knurl area) is made both on the body of the Watson Screwcap Tube and on each of well of the Piccolo Rack so that the Screwcap Tube is completely fixed. Thus, the cap of the Screwcap Tube can be easily opened and closed with one hand (much easier to open by using Screwcap Opener<#175-15C>).

Material: Polycarbonate

Cat. No.	Item	Dimensions	Unit	Case
1522-2451	Piccolo Rack for 24 Screwcap Tubes, Clear	118 × 82 × 52 mm	10 pieces	5 units
1522-2453	Piccolo Rack for 24 Microtubes, Clear	118 × 82 × 52 mm	10 pieces	5 units

Freezer Storage Racks



1577-550

- Excellent storage racks in terms of its characteristics, heat-proof and crash-proof.
- Conventional frozen preservation tube (to 2.0 mL) can be applied to and the plate can be piled up.
- Number and alphabet for identification are marked to the partition.
- Possible to store in the liquid nitrogen (vapor phase) and in the deep freezing condition.

Material: Polycarbonate (rack)  
Polypropylene (plate)



Cat. No.	Item	color	Dimensions	Unit	Case
1577-550B	Freezer Storage Rack for 50 Microtubes	Blue	146 × 76 × 55 mm	10 pieces	10 units
1577-550G	Freezer Storage Rack for 50 Microtubes	Green	146 × 76 × 55 mm	10 pieces	10 units
1577-550R	Freezer Storage Rack for 50 Microtubes	Red	146 × 76 × 55 mm	10 pieces	10 units
1577-550W	Freezer Storage Rack for 50 Microtubes	White	146 × 76 × 55 mm	10 pieces	10 units
1577-550Y	Freezer Storage Rack for 50 Microtubes	Yellow	146 × 76 × 55 mm	10 pieces	10 units
1577-850B	Freezer Storage Rack for 81 Microtubes	Blue	131 × 131 × 52 mm	10 pieces	10 units
1577-850G	Freezer Storage Rack for 81 Microtubes	Green	131 × 131 × 52 mm	10 pieces	10 units
1577-850R	Freezer Storage Rack for 81 Microtubes	Red	131 × 131 × 52 mm	10 pieces	10 units
1577-850W	Freezer Storage Rack for 81 Microtubes	White	131 × 131 × 52 mm	10 pieces	10 units
1577-850Y	Freezer Storage Rack for 81 Microtubes	Yellow	131 × 131 × 52 mm	10 pieces	10 units

Note:(\*)Repeated autoclaving may cause rack damaged.

## Cryotube Boxes

- The cryotube box for 40 tubes.

\* NOT Autoclavable



Material: Styrene resin

Cat. No.	Item	Dimensions	Unit	Case
1568-100	Cryotube Box for 40 tubes, 1.0 mL	125 × 83 × 42 mm	10 pieces	5 units
1568-200	Cryotube Box for 40 tubes, 2.0 mL	125 × 83 × 54 mm	10 pieces	5 units

## PCR Tube Racks

- PCR single tubes (0.2 mL) × 96 pieces can be stored. Frozen storage at -80°C. Autoclavable.

- The black-colored shading type is recommended for the storage of the light sensitive sample.



\*:see below



1512-952  
Material: Polycarbonate (rack)  
Polypropylene (plate)



1512-952BL  
Material: Polypropylene

Cat. No.	Item	Dimensions	Unit	Case
1512-952	PCR Tube Rack (PC) for 96 tubes	118 × 82 × 40 mm	10 pieces	10 units
1512-952BL	PCR Tube Rack (PP) for 96 tubes, Black for Shading	118 × 82 × 40 mm	10 pieces	10 units

Note:(\*)Repeated autoclaving may cause rack damaged.



**Centrifuge Tube Floating Rack “Jellyfish”**

**A centrifuge rack that can also be used as a floater.**



- Dedicated rack for four 15 mL / 50 mL centrifuge tubes.
- The rack can be used as a floater by changing the direction of the tray (See below).
- Made of polypropylene and can be autoclaved.



Material: Polypropylene


**Dual-use**

**As a rack**



Attach a tray to the bottom of the rack as a receptacle.

**As a floater**



If attaching a tray onto the rack body (after putting centrifuge tubes), the rack may be used as floater.

Cat. No.	Item	Unit
1541-809	Centrifuge Tube Floater Rack “Jellyfish”	1 piece

**Slide Box**



Material: Polypropylene

- 10 pieces of slides or counting plates can be compactly set.
- (Applicable size: 76 x 26 mm)

Note: Slides are not included.

Cat. No.	Item	Case
158-102BL	Slide Box for Cell Counting Plate / Microscope Slide, Black	30 pieces

Centrifuge Tube Floating Rack / Slide Box