

dermlite DL200



The DermLite DL200, in its beautifully ergonomic, easy-to-disinfect body, integrates DermLite's patented polarization system for the clearest views deep into the lesion, our curated color spectrum that highlights critical lesion structures, and super-fast changes between contact & non-contact dermoscopy. Features such as a removable glass faceplate for peripheral-contact dermoscopy, IceCap®

DL200HR

DermLite DL200 HR is permanently cross-polarized, with two brightness levels, making it the brightest polarized DL200.

What's unique

- + Always cross-polarized
- + Two brightness levels
- + Brightest DL200

You might consider

- DermLite DL200 Hybrid
- DermLite DL3N
- DermLite DL4

Polarization Indicator
Know your modality
at all times

disposable caps to comply with strict hygiene standards, hours of run time on a single charge and integrated Magneti-Connect® smartphone connectivity, and your choice of either exclusively polarized light in the *HR*, or both polarized & non-polarized modalities in the *Hybrid* model, combine to make the lightweight DL200 the most elegant & capable compact dermatoscope at your fingertips.

DL200 benefits

- + Compact, highly ergonomic design
- + Go between contact & non-contact with the flick of your finger
- + Same proven optics as DL2
- + Extremely long battery life
- + Mobile device-compatible
- + Four-level charge indicator
- + Faceplate with 10 mm reticle

Better devices may offer

Better optics, more precise focusing, more lighting options

DL200Hybrid

DL200 Hybrid combines polarized & non-polarized lighting. Visualize pigmented and vascular structures with or without immersion fluids and see the benefits of polarized and non-polarized light.

What's unique

- + Toggle between polarized & non-polarized
- + Available in silver & gold

You might consider

- DermLite DL3N
- DermLite DL4



IceCap®
Disposable snap-on
infection control



Removable Faceplate
Improved visualization of vascular
structures without contact

