



iPad



Tablets



DSLR



Compact Camera



Android Phones



iPhone



- + Polarized & non-polarized
- + Contact & non-contact dermoscopy
- + Superb 10x lens, large focal range
- + Fine-tuned blue frequency range of white light for improved rendition of vascular structures
- + PigmentBoost™ for pigmented lesions
- + Fully enclosed focusing mechanism
- + Camera/iPhone/iPad/Galaxy compatible
- + Faceplate with 10 mm markings
- + Convenient desktop charger

MagnetiConnect®

The quickest way to use DermLite with your smartphone or iPad. Sold separately.



DermLite DL3N and Nikon 1
Beautiful images in a compact package

dermlite DL3N

DermLite DL3N is the most capable third-generation DermLite dermatoscope. With fine-tuned illumination of amelanotic lesions plus PigmentBoost™ for optimized pigmented lesion imaging, it has been the tool of choice of the world's leading dermatologists.

Precision-engineered and crafted from recyclable aluminum, it features a 25 mm four-element lens, 28 high-powered LEDs and PigmentBoost™ illumination designed to see pigmented lesions better than ever before. For optimum ease of use, a tap on its power button instantly toggles between polarization modes, while a secondary button activates PigmentBoost.



Proven excellence.

The convenient DL3N retractable faceplate spacer design and focusing mechanism enables you to quickly switch between non-contact & contact modalities and lets you precisely focus your image while in skin contact. The spacer comes with a glass faceplate with 10 mm markings and, due to its fully enclosed design, contaminants as well as undesirable ambient lighting no longer affect your image. To reduce reflections and greatly improve the contrast of the image, an optional eyepiece is included. With the DermLite DL3N, it's easy to create your very own digital dermoscopy system by attaching one of many compatible cameras or camera-equipped mobile devices. Simply attach your digital



camera, smartphone or iPad, extend the spacer, and start taking sharp, glare-free images quickly and easily. This state-of-the-art device will let you know when it's time to recharge its long-lasting lithium battery: Once the charge level drops below 25%, the "3" on the top of device changes from green to orange. DermLite DL3N packs the unique advantages of both polarized and non-polarized dermoscopy with or without immersion fluids, improved vascular imaging, superior pigment visualization with PigmentBoost, great optics, solid durability, refined ergonomics and amazing versatility into a beautifully elegant design.