



Magnetically quick & simple  
DermLite DL4 can be magnetically attached to virtually any phone or tablet via a choice of DermLite adapters

- + compact all-aluminum design
- + patented polarization system for unsurpassed visual clarity
- + polarized & non-polarized
- + IceCap® infection control system
- + large 30 mm 10x lens system
- + charge from any USB port
- + PigmentBoost® Plus (not in DL4W)
- + contact & non-contact dermoscopy
- + magnetically attach smartphones & tablets (requires adapter)
- + 4 charge indicator levels
- + 2 brightness levels
- + snap-on faceplate
- + protective leather pouch included

Small Area Contact Plate  
Sold separately



## The most refined vision in dermoscopy.

The 4th-generation DermLite dermatoscope has been designed from the ground up to be more refined in every detail. An amazingly smooth, ergonomically optimized all-aluminum design that is as easy to disinfect as it comfortably molds to your hand, the ultra-compact DermLite DL4 puts phenomenal versatility, the largest optics, and DermLite's signature color spectrum that has been optimized to highlight the lesion structures most critical for early detection – right at your fingertips. DL4's patented polarization system provides a vision of unsurpassed clarity – deep into the lesion. Its glass faceplate (with 10 mm reticle) simply snaps on and your smart devices (smartphones & tablets) attach magnetically, enabling you to go between taking clinical and dermoscopic photos quicker than ever before.

Hygiene and infection control have become increasingly important aspects of providing care, and DL4's IceCap® has you covered: Disposable snap-on caps attach over the glass faceplate and may greatly reduce the risk of cross-contamination. A newly developed DL4 retractable glass faceplate design allows you to go between non-contact & contact modalities smoother than ever and lets you precisely focus your image while in skin contact. DermLite DL4's high-powered LEDs and rechargeable battery produce light in an extremely energy-efficient way. And, in order to know when you run out of battery,

DL4 now has a more granular four-level charge indicator. If you need to recharge, you can do so from any USB port using the standard Micro USB to USB cable included with each DL4.

**Versatile DL4 vs. Ultra-Bright DL4W**  
While model no. DL4W produces more light from its exclusively white LEDs, the top-of-the-line model no. DL4 provides enhanced imaging of deeper pigmentation with two levels of color spectrum control (PigmentBoost Plus).