

SYRIS v900L

Enjoy the World's Best Polarized Light, Hands-Free.

The Syris® v900L visualization system combines a proprietary superior white-light LED with hybrid polarization to enhance your vision by removing surface glare from the patient's skin, promoting improved viewing of sub-surface features.

The v900L is ideally suited for:

- Vascular laser treatments
- Surgical excisions & biopsies
- Sclerotherapy
- Hair-removal/grafting/transplantation
- Botox or filler treatments
- Tattoo removal




Without Syris



With Syris



Without Syris



With Syris

Ergonomic Headgear

Padded, dual-adjustment design for a most comfortable fit



Adjustable Polarization

Highlight surface and sub-surface features

Visi-Shield™

Unobstructed field of view, 2.0x and 2.5x magnification



Waist-Mounted Battery Unit

Dual brightness control in waist-mounted battery unit for up to 4 hours of use

Charger Included.

Convenient charging base allows you to charge a 2nd battery (sold separately) for nearly uninterrupted use.

