



WikiPouch®
by Infection Prevention Products 

CLEAN – SAFE – SECURE

Manufacturer Recommended Usage for the WikiPouch® formerly IP-Pouch

WikiPouch® Application Protocol effective December 2021

1. Write-in date when the WikiPouch® is installed. Use a permanent marker include Day/Month/Year.
2. The WikiPouch® is to be secured in close proximity to the patient for easy access.
3. The WikiPouch® can be adhered to oxygen concentrators, wheelchairs and bedside tables etc. Attaches 3 ways: Peel & Stick Tape, Drawstring, or Velcro Straps.
4. To protect the patients, staff & the facility from cross contamination, always store Nasal Cannulas & HHN's inside the WikiPouch® when not in use.
5. **The WikiPouch® is to be replaced once per month per patient or PRN.** (Replace when visibly soiled or if airflow is impeded).

Call if you have any questions:

Infection Prevention Products Inc.
DBA WikiPouch®
1124 Almond Vista CT.
Chico, CA 95973
www.IPPmfg.com
1-888-587-7775



A variety of WikiPouch® products attach in a range of useful applications

- Peel & Stick
- Drawstrings
- Velcro Straps