

Fact Sheet

Powered Mobility

Freerider Luggie Travel Power Chair

The Luggie Travel Chair folds with two simple actions into a package that can be easily stored when not in use. Even with the seat fully folded to the floor it can still be used as a traditional electric wheelchair, so very young people, some amputees, and others who would benefit from a very low seat can use this powerchair. It is suitable for indoor and gently sloping outdoor use.

* Portable, foldable design

Runner wheels to assist with moving when folded

Superior technology and innovative engineering

* IATA approved Li-ion batteries 10.5AH

3.1cm floor clearance

- * Max speed 6.4kph
- * 15-20km range
- * 5 amp charger
- Weight capacity 145kg



Specifications

Overall height	
Overall depth	95
Overall width	61
Seat depth	40
Seat width	43

Measurements in cm, weight in kg



Product Code: FRLCHAIR