

Medi Wool Overlay

The Australia Medical Sheepskin is used for health care, medical and institutional uses, meeting the Australian Medical Sheepskin Standard: AS4480.1 1998 = HiTemp UR Medical Sheepskin. By using a Medi Wool overlay it minimizes the incidence, severity and duration of pressure ulcers (bedsores) therefore decreasing friction, pressure and humidity at the points of contact between the patient and the bed or chair.

- * Offers excellent air redistribution properties of a wool product
- * Provide the ideal microclimate for a good sleep
- * Offers practical comfort absorbing over 36% of its dry weight in moisture
- * Looking after temperature and perspiration regulation
- * Keeps you warm during colder seasons and cool during summer
- * Washable at temperatures up to 80° C



Product Code: SS3001

Specifications

Size XL
Dimensions (mm) L1030 x W700