



Pressure Care



Vinyl Covered Gel Wheelchair Cushion

The Aidapt Gel cushion has a wipe-clean, vinyl cover and foam inner to provide optimum comfort to the user. The Gel layer keeps the cushion cool when using for long periods of time. The cushion can provide extra support and comfort in cars and homes.

- * Ideal for use in the home, car and on wheelchairs
- * Gel layer to keep cool whilst in use
- * Easy-care, removable cover is wipe clean vinyl
- * Product Dimensions 75x450x400mm
- * Net weight 1.54kg



Specifications

Product Code: VA126SC 45x40x7.5cm



Ph 1300 622 633 mobility@mountiescare.com.au www.mountiescaremobility.com.au