



Fisherbrand™ Non-sterile Transfer Pipettes, Transparent LDPE

edge»

Prevent cross contamination. Fisherbrand™ Non-sterile Transfer Pipets, Transparent LDPE are made of shatterproof, flexible plastic and are designed to replace glass Pasteur pipets. Graduated and non-graduated pipets available in a range of sizes and styles, including jumbo, thin stem and extra long.

Brand: Fisherbrand™

The Fisherbrand
Guarantee

Code: A0

Additional Details: Weight : 0.81000kg

Description

- Built-in bulb eliminates the need for plugging, packaging and autoclaving
- Flexible stem bends to sample from narrow or small tubes
- Low affinity surface reduces loss of cells and valuable proteins due to binding
- Soft plastic pipets made of non-toxic low-density polyethylene (LDPE)
- Non-sterile, transparent and fully recyclable

Specifications

Drops/mL	25	Graduated	Yes
Graduations	0.5 mL	Quantity	500 Pack
Tip Style	Standard	Volume (Metric) Dispense	40 µL drop
Length (Metric)	155 mm	Material	LDPE
Type	Transfer Pipet	Volume (Metric)	3 mL