

# Fisherbrand™ Disposable Plastic Syringes



These sterile single-use, disposable syringes are individually wrapped in a paper/plastic blister package. They are available in a wide range of capacities, have a concentric tip, and offer either Luer-lock or Luer-slip tips.

- Designed for small-volume filtration and aspiration of liquids
- Suitable for use with syringe filters
- Three-part construction includes a polypropylene barrel and piston and a rubber gasket
- For research purposes only

Cat. No	Capacity	Luer fitting	Pack qty
15889142	1mL	Luer-slip	100
15899142	3mL	Luer-slip	100
15809152	5mL	Luer-slip	100
15819152	10mL	Luer-slip	100
15829152	20mL	Luer-slip	50
15839152	60mL	Luer-slip	50
15849152	1mL	Luer-lock	100
15859152	3mL	Luer-lock	100
15869152	5mL	Luer-lock	100
15879152	10mL	Luer-lock	100
15889152	20mL	Luer-lock	50
15899152	60mL	Luer-lock	50

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20  
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478  
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00  
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303  
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500