


Apical Shaper[®]

NEXT-GENERATION INSTRUMENTATION
FOR APICAL SHAPING WITH 3% TAPER



AVALON
BIOMED.



Torque NCM 2-4 RPM: 500

Apical Shaper®

NEXT-GENERATION INSTRUMENTATION
FOR APICAL SHAPING WITH 3% TAPER



Apical Shaper's unique design gives you the best of both worlds. The parallelogram cross section delivers efficient cutting, while the .03 taper allows you to finesse the apical 1-2mm, which we know is the most important part of the root canal.

RYAN WALSH, DDS, MS
KELLER, TX

4 PACK ApicalShaper - \$28.⁰⁰ | 6 PACK ApicalShaper - \$39.⁰⁰

ApicalShaper 21mm

ZAS2103003 E30 21 mm (6 units) sterile
ZAS2103503 E35 21 mm (6 units) sterile
ZAS2104003 E40 21 mm (6 units) sterile
ZAS2105003 E50 21 mm (6 units) sterile
ZAS21ASST1 E30-E50 21 mm (4 units) sterile

ApicalShaper 25 mm

ZAS2103003 E30 25 mm (6 units) sterile
ZAS2103503 E35 25 mm (6 units) sterile
ZAS2104003 E40 25 mm (6 units) sterile
ZAS2105003 E50 25 mm (6 units) sterile
ZAS21ASST1 E30-E50 25 mm (4 units) sterile

ApicalShaper 31 mm

ZAS2103003 E30 31 mm (6 units) sterile
ZAS2103503 E35 31 mm (6 units) sterile
ZAS2104003 E40 31 mm (6 units) sterile
ZAS2105003 E50 31 mm (6 units) sterile
ZAS21ASST1 E30-E50 31 mm (4 units) sterile