

SlimShaper®

THE FIRST MINIMALLY
INVASIVE THREE-ALLOY SYSTEM



AVALON
BIOMED.

	MFD	TAPER	Ø
ZS1	0.80	2-6%	0.15
ZS2	0.80	4-3%	0.20
ZS3	0.84	4-3%	0.25

GOLD - Maximum cutting capacity and torsional resistance.

PINK - The perfect balance between cutting capacity and flexibility

BLUE - Extraordinary flexibility and greater cyclic fatigue resistance

Torque NCM 2-4 RPM: 500

SlimShaper®

THE FIRST MINIMALLY
INVASIVE THREE-ALLOY SYSTEM



My cone fit radiographs are consistent and rarely need any refinement. I cannot recommend SlimShaper enough. It's one of the best MIE files in its class.

SAMEER JAIN, DDS, MS, MSD
WOODLAND HILLS, CA

6 PACK SlimShaper - \$42.00

SlimShaper® 21 mm

ZSS2100001 ZS1 21 mm (6 units) sterile
ZSS2100002 ZS2 21 mm (6 units) sterile
ZSS2100003 ZS3 21 mm (6 units) sterile
ZSS21ASST6 ZS1-ZS3 21 mm (6 units) sterile

SlimShaper® 25 mm

ZSS2500001 ZS1 25 mm (6 units) sterile
ZSS2500002 ZS2 25 mm (6 units) sterile
ZSS2500003 ZS3 25 mm (6 units) sterile
ZSS25ASST6 ZS1-ZS3 25 mm (6 units) sterile

SlimShaper® 31 mm

ZSS3100001 ZS1 31 mm (6 units) sterile
ZSS3100002 ZS2 31 mm (6 units) sterile
ZSS3100003 ZS3 31 mm (6 units) sterile
ZSS31ASST6 ZS1-ZS3 31 mm (6 units) sterile

SlimShaper® Gutta-percha

ZSSGP00ZS2 Gutta Percha ZS2 (60 units)
ZSSGP00ZS3 Gutta Percha ZS3 (60 units)

SlimShaper® Paper Points

ZSSPP00ZS2 Paper Points ZS2 (60 units)
ZSSPP00ZS3 Paper Points ZS3 (60 units)