



Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 33,000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 5 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.

In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



www.nipro-group.com/en-en/our-company/our-locations
Please contact your local representative for more information.




Nipro Australia Pty Ltd : Suite 2.02, Level 2, 657 Pacific Hwy, St Leonards NSW 2065 | T: +61 2 9437 9000 or 1800 451 737 | info@niproaustralia.com.au | niproaustralia.com.au

Nipro Medical Europe : European Headquarters, Blokhuisstraat 42, 2800 Mechelen, Belgium | T: +32 (0)15 263 500 | www.nipro-group.com

Nipro Corporation Japan : 3-9-3 Honjo-nishi, Kita-ku, Osaka 531-8510, JAPAN | T: +81 6 6372 2331 | www.nipro.co.jp

Products are continuously under review in the light of technical advancement; the actual specifications may therefore be subject to improvement or modification without notice. Some features and/or accessories designed and described could be applied in different markets. The information provided does not constitute an offer to sell, and is provided for merely informational purposes. Images published in the present document are shown exclusively for demo purposes and might be adapted for specific needs. All rights reserved.

MM-1696 - Bro-Acid_Concentrates - EN - 18.Feb.21



RENAL CARE

ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS



ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS

Ready-to-use liquid acid concentrates intended for use with bicarbonate concentrate.

Acid concentrate is an essential component in the preparation of dialysis fluid, whereby the hemodialysis machine mixes liquid acid concentrate with bicarbonate concentrate and purified water.

The standard mixing ratio of Nipro acid concentrates is 1:44.

PORTFOLIO

	Na ⁺ mmol/l	Ca ⁺⁺ mmol/l	Mg ⁺⁺ mmol/l	K ⁺ mmol/l	Cl ⁻ mmol/l	CH ₃ COO ⁻ mmol/l	HCO ₃ ⁻ mmol/l	Glucose g/l	Osmolarity mOsm/l
306-A (A+B)	139	1,25	0,48	1,50	103,14	2,85	37,96	1,50	294,50
307-A (A+B)	139	1,75	0,48	1,50	104,14	2,85	37,96	1,50	296,00
308-A (A+B)	139	1,50	0,48	1,50	103,64	2,85	37,96	1,00	292,48
310-A (A+B)	139	1,50	0,50	1,00	103,20	2,85	37,96	1,00	291,56
313-A (A+B)	140	1,50	0,50	1,50	108,50	3,00	34,00	1,00	294,55
314-A (A+B)	140	1,25	0,50	1,50	108,00	3,00	34,00	1,00	293,80
316-A (A+B)	140	0,63	0,25	1,50	102,44	4,00	36,81	1,00	291,18
318-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,50	299,08
319-A (A+B)	139	1,75	0,48	3,00	105,64	2,85	37,96	1,50	299,00
320-A (A+B)	139	1,25	0,50	3,00	104,70	2,85	37,96	1,50	297,59
321-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,00	296,30
323-A (A+B)	140	1,50	0,50	2,00	109,00	3,00	34,00	1,00	295,55
324-A (A+B)	140	1,25	0,50	2,00	108,50	3,00	34,00	1,00	294,80
326-A (A+B)	139	1,50	0,48	2,50	104,64	2,85	37,96	1,50	297,25
327-A (A+B)	140	1,50	0,50	3,00	110,00	3,00	34,00	1,00	297,55
330-A (A+B)	140	1,75	0,50	3,00	110,50	3,00	34,00	1,00	298,30
331-A (A+B)	140	1,25	0,50	3,00	109,50	3,00	34,00	1,00	296,80
332-A (A+B)	140	1,50	0,50	4,00	111,00	3,00	34,00	1,00	299,55

PACKAGING DETAILS

	Volume (L)	Reference code	Units per box	Units per pallet
Bags	3,8	BSOL45-*A4LBX	3	144
	4,7	BSOL45-*A5LBX	2	128
Stackable canister(s)	3,9	SOL45-*A4L	N/A	132
	5	SOL45-*A5L	N/A	132
	10,5	SOL45-*A11L	N/A	60
Canister(s) in a box	3,9	SOL45-*A4LBX	2	128
	5	SOL45-*A5LBX	2	128
Octabin - Bottom connector	500	CSOL45-*A500L	N/A	1
	300	CSOL45-*A300LT	N/A	1
Octabin - Top connector	500	CSOL45-*A500LT	N/A	1
	600	CSOL45-*A600LT	N/A	1
Big box stackable canisters	3,9	SOL45-*A4LBB	96	96
	5	SOL45-*A5LBB	96	96

* acid concentrate formula

Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Please contact your Nipro representative for product availability and information.

SHELF LIFE: 2 years

