


# DIAMOND KNIVES

**NATURAL DIAMONDS  
UP TO 3000 INCISIONS**

## CATARACT SURGERY

### SIDE PORT KNIVES



Universal Straight Titanium Handle, 110 mm


6-10/6-050

45° Single Edge  
1.00 mm



6-10/6-051

30° Single Edge  
1.00 mm



**MOST  
POPULAR**

6-10/6-052

45° Double Edge  
1.00 mm



**MOST  
POPULAR**

6-10/6-053

Trifacet  
1.00 mm



6-10/6-0531

20° Trifacet  
1.00 mm



6-10/6-070

60° Lancet  
1.00 mm



6-10/6-0701

60° Lancet  
0.80 mm



6-10/6-076

Multi-Incision  
1.00 mm



### PHACO KNIVES



Universal Angled Titanium Handle, 130 mm

#### Clear Cornea Blades

6-20/6-075 – 2.20 mm **MOST POPULAR**  
 6-20/6-071 – 2.50 mm  
 6-20/6-077 – 2.75 mm  
 6-20/6-080 – 2.65 mm  
 6-20/6-072 – 2.80 mm **MOST POPULAR**  
 6-20/6-073 – 3.00 mm  
 6-20/6-074 – 3.20 mm



#### Trapezoid Self-Diving Blades

6-20/6-143 – 1.50/1.80 mm  
 6-20/6-144 – 1.80/2.00 mm  
 6-20/6-145 – 1.80/2.20 mm  
 6-20/6-107 – 2.00/2.30 mm **MOST POPULAR**  
 6-20/6-100 – 1.90/2.50 mm  
 6-20/6-104 – 2.30/2.80 mm **MOST POPULAR**  
 6-20/6-105 – 2.60/3.00 mm  
 6-20/6-101 – 2.70/3.20 mm  
 6-20/6-102 – 3.00/3.50 mm



#### Zaldivar Blade

6-20/6-0551 – 0.55/1.00 mm  
 Designed for ICL™ implantation




#### Facet-Free Crescent Blades

6-20/6-091 – 1.50 mm  
 6-20/6-092 – 2.00 mm  
 Can be used for glaucoma surgery



## LRI SURGERY

### UNIVERSAL THREE-STEP LRI KNIFE



**6-322/6-0531**  
 20° Trifacet Blade  
 0.20–1.00 mm  
 Single footplate for  
 better visualization




Titanium Handle with Calibration Mechanism, 130 mm

- Precise, easy and stable calibration mechanism with 3 depth settings of 500, 550 and 600 microns
- Full extension for side port incision

### PRE-SET LRI KNIVES

20° Trifacet Blade  
 0.20–1.00 mm



**6-500/6-0531**  
 Preset depth 500 microns

Standard Handle, 110 mm

**6-600/6-0531**  
 Preset depth 600 microns

**6-500S/6-0531**  
 Preset depth 500 microns


**6-600S/6-0531** **MOST POPULAR**  
 Preset depth 600 microns



Short Handle, 78 mm

## AK & ICRS PROCEDURES

### CALIBRATED DIAMOND KNIVES




Titanium Micrometer Handle, 105 mm


**MOST POPULAR**

**6-00/6-020**

- Performs a 1.00 mm radial incision at 70-80% of corneal thickness
- Designed for the implantation of intracorneal ring segments




**6-00/6-030**  
 Tapered blade



## CLEANING PACK

### DIAMOND KNIFE CLEANING PACK

The pack consists of three solutions that are used to clean the blade of residual debris prior to sterilization.



21-602-1