


Luxurious adjustable bed collection for an enhanced sleeping experience





ComfiMotion collection

We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.

Create Your ComfiMotion





Bed Function

Headboard Style








Accessories





Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Auto-Regression in the Backrest: Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution.

Clinical Positioning: One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

Upholstery



Portobello Grey

Leg Raise: Promotes circulation, relieving tired, swollen/ fatigued or painful limbs.

Memory Function: Allows the user to set and return to their preferred positions at any time with one click of a button.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.

Lockable Castors: Eight heavy duty dual castors.

*Carers should seek advice from a trained Allied Health professional for



▲ Companion Sleeping See Configurations sec on for more companion options

Sizes



King Single



Features

- risk, improving user safety.
- overbed tables.
- during power outages.
- visible at all times.
- and back relief.

Supporting the modern, independent lifestyle.

The Activ Care bed combines the best of in-home care functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

• Ultra Low Height: Minimum platform height of 195mm reduces falls

• Customisable Safe Entry/Exit Height: Programmable single button press ensures perfect bed height to maximise independence.

• Underbed Clearance: Provides clear access for patient lifters and

• Battery Backup: Provides confidence and maintains independence

• Night Assist: Luminous underbed halo for safe entry and exit.

• Illuminated Handset: Large button interface is intuitive and highly

• Massage: Relaxing vibration modes for personalised neck, shoulder




2030 x 1070 x 195 – 700mm



Visit comfimotion.com.au for Limited Warranty T&C's

Care 🖄 Adjustable Bed



Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Clinical Positioning: Raises back and legs to alleviate cardiac exertion, improving resting comfort.

Leg Raise: Promotes circulation, relieving tired, swollen/ fatigued or painful limbs.

Upholstery

Portobello Grey

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.

Lockable Castors: Four heavy duty castors.

*Carers should seek advice from a trained Allied Health professional for

STANDARD

350 – 650mm with

Standard Castors

Variable Bed Height Options



280 – 580mm with Low Profile Bed Legs



ACCESSIBILITY 400 - 700mm with Large Castors (optional)





▲ Companion Sleeping

Sizes





1. With Low Profile Bed Legs.

6

For comfort you deserve.

The Care adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.

- Night Assist: Luminous underbed halo for safe entry and exit
- Enclosed Frame Design: Ensures optimal positioning when getting
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo
- Box-on-Box Design: Facilitates easier entry and exit and allows true

2030 x 1070 x 280 – 580 mm¹



Visit comfimotion.com.au for Limited Warranty T&C's



Pocket Spring

The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality.




Transitional Layered Pocket Springs arranged in weight activated support zones.




Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



Luxuriously Soft Pillowtop Layer provides a plush immersive feeling complimented by the responsive pocket springs.




Our ComfiMotion Pocket Spring and Memory Foam mattresses are designed with the flexibility to follow the contouring profile of an adjustable bed.



Classic Pocket Spring Mattress

Pocket Spring



2030 x 1070 x 240mm





Memory Foam



King Single – Classic 2030 x 1070 x 200mm **King Single – Incontinence** 2030 x 1070 x 200mm

Memory Foam

The ComfiMotion Memory Foam mattress contours to the user's body, providing personalised comfort and transitional pressure relief. Integrated sidewalls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit.



Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.



Superior Breathability: Up to 8x more breathable than conventional memory foams.



Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.



Eco-friendly GECA Foam is produced using cleaner chemicals and production technology.



▲ Classic Memory Foam Mattress




Visit comfimotion.com.au for Limited Warranty T&C's

Configurations for Companions

Whichever ComfiMotion package you select, rest assured there is a companion solution available to suit those wishing to sleep alongside any ComfiMotion Adjustable bed.







Dual Adjustable Beds

Designed to sit flush when positioned together to enable close monitoring and enhanced care. This configuration offers the flexibility of independent adjustability, without compromising on the sense of inclusivity and intimacy between partners.

Split King Upgrade Kit

Ensures simultaneous Hi–Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.

Split King Upgrade Kit

Fixed Companion Bed - (LS) 2030 x 920 x 370mm

Fixed Companion Bed

Matching single bed without any

adjustability designed to sit flush

when positioned together. It facilitates


a side-by-side sleeping arrangement

at a fixed height of 370mm, without

compromising on the users care and

assistance requirements.


Headboards




BEA830100 - Long Single BEA830110 - King Single BEA830116 - King

Allure BEA830120 - Long Single BEA830130 - King Single BEA830136 - King

Accessories



Collection Summary



Pillows







Plush

Large Accessibility Castors MWT451080


(suits Care Bed only, 4 x required)

Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access

Battery Backup BEA830220 - Care Bed *Standard on Activ Care Bed

Provides peace of mind in the event of a power outage.





Simplified Handset BEA830214 - Activ Care Bed

For clients with reduced cognitive ability.

Split King Upgrade Kit BEA830230

Ensures simultaneous Hi–Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.

Headboard Headboard		Allure Headboard – (LS)
		Allure Headboard – (KS)
		Allure Headboard – (K)
		Chesterfield Headboard – (LS)
		Chesterfield Headboard – (KS)
	BEA830116	Chesterfield Headboard – (K)
	BEA830200	Crook Handle Bed Stick – Care Bed
		Crook Handle Bed Stick – Activ Care Bed
		Bed Rail – Care Bed
		Bed Rail – Activ Care Bed
		Large Accessibility Castors – Care Bed
		Battery Backup – Care Bed
		Simplified Handset – Activ Care Bed
		Split King Upgrade Kit – Activ Care Bed
	Key: LS - Long	g Single KS – King Single K – King

Clinical Considerations

Accessibility and Inclusivity: Unexposed bed surround with padded upholstery removes gaps between beds in co–sleeping arrangements and promotes safe entry and exit.

User Independence: Height adjustability improves user mobility and assists carers with manual handling.

Flexible Positioning: Caters to various user requirements such as respiratory conditions, snoring, sleep apnoea, pressure redistribution, joint or muscle pain, lymphoedema and swelling.

Vision and Hand Function: Handset with a clear and intuitive button display, assists users with reduced vision or hand function.

Falls Prevention: Integrated bed lighting and a minimum platform height of 195mm (Activ Care) and 280mm (Care) helps to reduce the likelihood and severity of falls in the home.

User Safety: Battery backup provides peace of mind in the event of a power outage (Standard on Activ Care, optional accessory on Care).



1111111111111

ASBR28.66 03032