

Enable Lifecare PTY LTD ABN 86 620 684 586 ACN 620 684 586 support@enablelifecare.com.au 1300 370 370

enablelifecare.com.au

CONFIGURA COMFORT ASSESSMENT FORM

Number: Email: Funding Body	/ :		A				
Chair Size	Part Number		Seat Depth	Seat Height	Seat Width		Weight Capacity
☐ SMALL	☐ CR5404 - Duratek ☐ CR5422 - Vinyl		406 mm (16") #1 457 mm (18") 508 mm (20")	388 mm (15") #4 406 mm (16") 431 mm (17") #2 457 mm (18") #3	457 mm (18")		130 kg
☐ MEDIUM	☐ CR5406 - Duratek☐ CR5424 - Vinyl		406 mm (16") #1 457 mm (18") 508 mm (20")	439 mm (17") #4 457 mm (18") 482 mm (19") #2 508 mm (20") #3	508 mm (20")		160 kg
☐ LARGE	☐ CR5408 - Duratek		406 mm (16") #1 457 mm (18") 508 mm (20")	439 mm (17") #4 457 mm (18") 482 mm (19") #2 508 mm (20") #3	559 mm (22")		160 kg
Kit 2" (50 mm) Kit 1"		ight Adjustment (25 mm)	#3 Height Adjustmer Kit 2" (50 mm) CA2422	nt 🗆	#4 Height Reduction Kit 1" (18 mm) CA2423		
PRESSON PI	essure			Post	ure		















Enable Lifecare PTY LTD
ABN 86 620 684 586 ACN 620 684 586
support@enablelifecare.com.au
1300 370 370

enablelifecare.com.au

Transfer Options

Description	Long Description	Part Number	Select
Drop-down Armside Left-hand side (as sitting in the chair)	Transfers; Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek CA2484 Small - Vinyl CA2492 Medium / Large - Duratek CA2488 Medium / Large - Vinyl CA2496	S
Drop-down Armside Right-hand side (as sitting in the chair)	Transfers; Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek CA2486 Small - Vinyl CA2494 Medium / Large - Duratek CA2490 Medium / Large - Vinyl CA2498	S



Extras



Cognitive



Support

Armrests	
Overhang on outside (std)	
Overhang on inside - reduces width between the armrests by 1.5" (38mm)	



Legrest Chan	nel	
Small CA2520		
Medium CA2521		
Large CA2522		9

Rise Setting			Legrest Setting		Handset Control	
Rise and Tilt (std) 23°	160kg SWL*		Standard		RHS (std)	
Medium Rise 14°	135kg SWL*		Standard + 1" (25 mm)		LHS	
Flat Rise 2°	110kg SWL*		Standard + 2" (50 mm)		Carer position	

^{*}MEDIUM & LARGE CHAIRS

Chair Assessment and Delivery Environmental Questionnaire

When carrying out an assessment please fill out the below form prior to ordering the chair.

Please give as much detail as possible.

ACCESS

Is there parking/access for a large van?	
Are there steps to the property?	
Is it easy access to the property?	
INTERNAL	
Are the doors and corridors wide enough for a chair to go through? Min. width for most chairs is 650mm.	
Does the chair have to go upstairs? If yes, please specify number of steps and if a stair lift is fitted.	
Is space being cleared for the chair?	
Is there a suitable power supply for the chair?	
INFECTION CONTROL	
Are there pets? Is the house smoke free?	

OTHER

Does the driver need to be aware of any other points?