

## REHABASSIST Slide Assist

	<p><b>Easy Slide Flat Slidesheet</b> - Polyester (<i>145cm x 200cm</i>) With the sliding surface inward, use 2 sheets (one of top of the other) or fold one sheet over to slide someone across a bed (from folded side to open side).</p> <p style="text-align: right;">SP145200FB</p>
	<p><b>Easy Slide Flat Slidesheet + Nylon handles</b> - Polyester (<i>90cm x 200cm</i>) Continuous handles around 3 sides of the sheet for easy grip. Use 2 together with sliding surface inward, or with a normal Easy slide flat sheet.</p> <p style="text-align: right;">S90200HO</p>
	<p><b>Positioning Slidesheet</b> - Open - Cushioned Polyester (<i>137cm x 46cm</i>) Use to slide someone into the back of their chair, or to facilitate moving them across surfaces using a transfer board.</p> <p style="text-align: right;">PS13746SO</p>
	<p><b>Gloves</b> - per pair - Polyester One carer can move someone in three steps - ankles, buttocks, shoulders (head on cushion). Two carers can slide someone up a bed if each has a glove. Can place under heels to reduce shear.</p> <p style="text-align: right;">G5120S</p>
	<p><b>One Way Slide</b> - Tube - Blue Seal Non-slip - (<i>37cm x 43cm</i>) Non-slip for use on slippery or woven surfaces. Position a person in the back of their seat and prevent slipping forward.</p> <p style="text-align: right;">O3743NT</p>
	<p><b>One Way Positioning Slide</b> - Open - Blue Seal Non-slip - (<i>130cm x 40cm</i>) on-slip long single layer. Can be folded with the handle placed through the back of an open-back chair or wheelchair and the carer steps on the handle to rather use their strong leg muscles to shift the person.</p> <p style="text-align: right;">O13040NO</p>
	<p><b>One Way Slide</b> - Open - Velour - (<i>37cm x 43cm</i>) Slightly padded velour one way slide. Help position person on pressure cushions or for positioning on a swivel cushion.</p> <p style="text-align: right;">O3743VO</p>
	<p><b>One Way Slide</b> - Tube - Velour - (<i>37cm x 43cm</i>) Velour cover with a light fill for added comfort. Only suitable for use on woven fabric surfaces so the fibres can grip.</p> <p style="text-align: right;">O3743VT</p>
	<p><b>One Way Positioning Slide</b> - Open - Velour - (<i>130cm x 40cm</i>) Non-slip long single layer for use on woven surfaces. Use folded with the handle through the back chair or wheelchair.</p> <p style="text-align: right;">O13040VO</p>