

FLOWFLEX™ Rapid SARS-CoV-2 antigen detection test

Metzcare reference: L031-11815

The Flowflex rapid test™ s SARS-CoV-2 is a rapid diagnostic test (TROD) for professional use only that detects COVID-19 by performing a single test on a patient*



- Easy to use, results in 15 minutes maximum
- Nasal or nasopharyngeal swab
- Test referenced by the ANSM and approved by the Ministry of Health:
<https://covid-19.sante.gouv.fr/tests>

Manufacturer/Laboratory	Acon Biotech
Declaration of CE conformity	Download
ISO 13485	Download
Flyer for display	Download
User manual FR	Download
Packaging	25 kits per box, 32 boxes per carton // 800 kits per carton 530*445*410mm, 13kg/carton
Box contents	25 test cassettes 25 sterile nasal or nasopharyngeal swabs 25 extraction buffer tubes 1* Positive and negative control swabs 1* Instruction manual 1* Sample collection guide
Clinical data	Clinical sensitivity: 97.60% Clinical specificity: 99.4%



Distributed by: Metze SAS
10 rue de Penthievre, 75008 Paris
commandes@metzecare.com



ACON Biotech (Hangzhou) Co., Ltd.