

Magnum-compatible biopsy needle



General Features

Manufacturer	Wellgo Medical Products GmbH
Country of origin	Germany
CE Mark	Medical Device – Class IIa European Directive 93/42/EEC
Sterility	STERILE (ISO 11737-2) ETO (ISO 10993-7)
Materials	Needle: S/S304 Protection: LDPE
Use	Disposable, single use
Shelf life	5 years
Advantages	Possibility of taking multiple samples from the same patient
	Echogenic needle tip for optimal ultrasound visibility
	Color coding available for easy identification of needle sizes
	Cannula marked in centimeters
	Compatible with W-Core II and Bard Magnum reusable guns

Pyrogen-free product: Yes

Latex-free: No

Reference table

REFERENCE	09cm	10cm	15cm	16cm	20cm	25cm	30cm
14G	021409W	021410W	021412W	021416W	021420W	-	-
16G	021609W	021610W	021615W	021616W	021620W	021625W	-
18G	021809W	021810W	021815W	021816W	021820W	021825W	021830W
20G	022009W	022010W	022015W	-	022020W	-	-

Packaging and wrapping

Detailed packaging	
Individual packaging	Master carton 20 pieces
Paper/plastic pouch	48*10*12cm 2 kg



Imported & Distributed by: Metze SAS

10 rue de Penthièvre, 75008 Paris

commandes@metzecare.com