

# Robotic Pipette Tips

---



**Metze**  
*care*

**PakGent**® 39/40

PakGent Bioscience offers a wide range of pipet tips for automated robotic workstations like, Tecan, Hamilton, Beckman etc.

All of our sterile pipette tips are manufactured out of the highest grade resins and are certified to be RNase, DNase, DNA, pyrogen, and ATP free. Our non-sterile tips are made from the same high quality resins and are certified RNase and DNase free as well.

Automated workstations actually require much tighter tolerances to the robotic tips than those of standard manual pipet tips. So we strongly suggest to purchase pre-sterilized products than non sterilized products which needs autoclaving, and autoclaving can warp the tips slightly.



**Robotic Tip Tecan®**

Fits Tecan Genesis Freedom®, Freedom EVO® and Miniprep with LiHa®

**50uL Clear Tip, sterile**




Part NO.	With Filter	Low Retention	Package Type	Pack
R 050	-	-	rack pack	96/rack 30racks/case
R 050 -L	-	Yes	rack pack	96/rack 30racks/case
R 050 F	Yes	-	rack pack	96/rack 30racks/case
R 050 F-L	Yes	Yes	rack pack	96/rack 30racks/case

**200uL Clear Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
R 200	-	-	rack pack	96/rack 30racks/case
R 200 -L	-	Yes	rack pack	96/rack 30racks/case
R 200 F	Yes	-	rack pack	96/rack 30racks/case
R 200 F-L	Yes	Yes	rack pack	96/rack 30racks/case

**1000uL Clear Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
R 1000	-	-	rack pack	96/rack 60racks/case
R 1000 -L	-	Yes	rack pack	96/rack 60racks/case
R 1000 F	Yes	-	rack pack	96/rack 60racks/case
R 1000 F-L	Yes	Yes	rack pack	96/rack 60racks/case

**50uL Conductive Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
R 050 X	-	-	rack pack	96/rack 30racks/case
R 050 X-L	-	Yes	rack pack	96/rack 30racks/case
R 050 FX	Yes	-	rack pack	96/rack 30racks/case
R 050 FX-L	Yes	Yes	rack pack	96/rack 30racks/case

**200uL Conductive Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
RT200TX	-	-	rack pack	96/rack, 30racks/case
RT200TX-L	-	Yes	rack pack	96/rack, 30racks/case
RT200TFX	Yes	-	rack pack	96/rack, 30racks/case
RT200TFX-L	Yes	Yes	rack pack	96/rack, 30racks/case

**1000uL Conductive Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
RT1000TX	-	-	rack pack	96/rack, 60racks/case
RT1000TX-L	-	Yes	rack pack	96/rack, 60racks/case
RT1000TFX	Yes	-	rack pack	96/rack, 60racks/case
RT1000TFX-L	Yes	Yes	rack pack	96/rack, 60racks/case

**Robotic Tip Hamilton®**

Fits Hamilton Star®, Starlet and Nimbus®

**50uL Clear Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
RT050H	-	-	rack pack	96/rack, 30racks/case
RT050H-L	-	Yes	rack pack	96/rack, 30racks/case
RT050HF	Yes	-	rack pack	96/rack, 30racks/case
RT050HF-L	Yes	Yes	rack pack	96/rack, 30racks/case

**300uL Clear Tip, sterile**



Part NO.	With Filter	Low Retention	Package Type	Pack
RT300H	-	-	rack pack	96/rack, 30racks/case
RT300H-L	-	Yes	rack pack	96/rack, 30racks/case
RT300HF	Yes	-	rack pack	96/rack, 30racks/case
RT300HF-L	Yes	Yes	rack pack	96/rack, 30racks/case

1000uL Clear Tip, sterile



Part NO.	With Filter	Low Retention	Package Type	Pack
T1000H	-	-	rack pack	96/rack, 60racks/case
RT1000H-L	-	Yes	rack pack	96/rack, 60racks/case
RT1000HF	Yes	-	rack pack	96/rack, 60racks/case
RT1000HF-L	Yes	Yes	rack pack	96/rack, 60racks/case

50uL Conductive Tip, sterile



Part NO.	With Filter	Low Retention	Package Type	Pack
RT050HX	-	-	rack pack	96/rack, 30racks/case
RT050HX-L	-	Yes	rack pack	96/rack, 30racks/case
RT050HFX	Yes	-	rack pack	96/rack, 30racks/case
RT050HFX-L	Yes	Yes	rack pack	96/rack, 30racks/case

300uL Conductive Tip, sterile



Part NO.	With Filter	Low Retention	Package Type	Pack
RT300HX	-	-	rack pack	96/rack, 30racks/case
RT300HX-L	-	Yes	rack pack	96/rack, 30racks/case
RT300HFX	Yes	-	rack pack	96/rack, 30racks/case
RT300HFX-L	Yes	Yes	rack pack	96/rack, 30racks/case

1000uL Conductive Tip, sterile



Part NO.	With Filter	Low Retention	Package Type	Pack
RT1000HX	-	-	rack pack	96/rack, 60racks/case
RT1000HX-L	-	Yes	rack pack	96/rack, 60racks/case
RT1000HFX	Yes	-	rack pack	96/rack, 60racks/case
RT1000HFX-L	Yes	Yes	rack pack	96/rack, 60racks/case

Robo c Tip Beckman®

Fi s Beckman FX/NX, Mult mek AP 6 and B omek3000

50uL Clear Tip, rack, sterile



Par NO.	W h Filter	Low Retent on	Package Type	Pack
T0 0B	-	-	rack pack	96/rack, 0racks/case
T0 0B-L	-	Yes	rack pack	96/rack, 0racks/case
T0 0BF	Yes	-	rack pack	96/rack, 0racks/case
T0 0BF-L	Yes	Yes	rack pack	96/rack, 0racks/case

250uL Clear Tip, rack, sterile



Par NO.	With Filter	Low Retention	Package Type	Pack
RT250B	-	-	rack pack	96/rack, 0racks/case
RT250B-L	-	Yes	rack pack	96/rack, 0racks/case
RT250BF	Yes	-	rack pack	96/rack, 0racks/case
RT250BF-L	Yes	Yes	rack pack	96/rack, 0racks/case

Robo c Tip Agilent®

Fi s Ag len ®, Ag len Bravo®

250uL Clear Tip, rack, sterile



Par NO.	With Filter	Low Retention	Package Type	Pack
T2 0A	-	-	rack pack	96/rack, 0racks/case
T2 0AF	Yes	-	rack pack	96/rack, 0racks/case

