



Professional  
**OEM/ODM**  
 Manufacturer  
 Health Care Hardware  
 Home Use Medical Devices



**Personal Health Care Hardware  
 Home Medical Electronics**

- Muscle Conditioning
- Pain Relief
- Massage & Relax
- Wrap Support



[www.leadingind.com](http://www.leadingind.com)

**Shenzhen Leading Technology Co.,Ltd**

TEL.: +86 755 2999 1449      Email: [info@leadingind.com](mailto:info@leadingind.com)

Web: [www.leadingind.com](http://www.leadingind.com)

Add.: 6F, Building C, No.5, Third Industrial Park, Jiuwei Community, Hangcheng Street, Bao'an District, Shenzhen, Guangdong China.

**Shenzhen Leading Technology Co.,Ltd**



# About Us



Established in 2014, Shenzhen Leading Technology Co., Ltd. is a professional OEM/ODM manufacturer integrating industrial design, R&D and manufacturing in the category of personal health care hardware and home medical electronics.

The company has always adhered to the development concept of "taking creative research and development as the core and closed-loop manufacturing as the foundation". It mainly focuses on technological innovation in the two fields of electric pulse (EMS&TENS) and Red Light Therapy (RLT), providing users with head-to-toe solutions and keeping to create high-quality as well as innovative health care hardware and home use medical devices.

The company consists of 1 R&D center, 5 production workshops, 6 production lines, more than 100 employees, and 6,000 square meter factory area. It has achieved one-stop closed-loop supply chain capabilities from mold manufacturing, plastic injection, hot pressing, fabric sewing, SMT to finished product assembly. Through excellent industrial design capabilities and rich experience in R&D as well as rapid response support of low-cost supply chain, we successfully incubated our own brand VONMIE in the Asia-Pacific market in 2017. Meanwhile, those capabilities greatly boosted sustainable development of some D2C brands held by our strategic customers in other markets.

Our products have various forms and solutions towards different parts of user body. Among them, the belt series has become a new business that has grown rapidly in recent years. The core products mainly include Heated Ab Toning Belt Stimulators, Red Light Pain Relief Belts, 2 in 1 Lower Back Pain Relief Brace Belt Stimulators, etc. So far, all of products cover more than 20 countries and regions around the world. Our overseas markets mainly include the United States, the European Union, Japan, Russia, South America and some Southeast Asian countries. Product certifications include but are not limited to FDA, FCC, CE, PSE, KC, REACH, ROHS, ANVISA, INMETRO, etc., which are fully compliant with the requirements of local laws and regulations in the above markets.

As a socially responsible manufacturer, our factory has obtained ISO9001 and ISO13485 quality management system and BSCI certification. We are committed the sustainable development plan to ensure our working conditions suitable for all employees and friendly to the environment.

All employees are committed to serving customers and increasing customer satisfaction through continuous improvement of product and service quality. We will adhere to the belief of "customers priority, employees first, and mutual development" and strive to work with our partners to realize self-worth, corporate value and social value.




# Contents

- 01 / EMS & TENS & RLT Introduction ..... 01-02
- 02 / Muscle Conditioning ..... 03-14
- 03 / Pain Relief ..... 15-24
- 04 / Massage & Relax ..... 25-28
- 05 / Wrap Support ..... 29-32


# EMS & TENS

EMS (Electrical Muscle Stimulation) and TENS (Transcutaneous Electrical Nerve Stimulation) involve applying electrodes to the skin near an affected muscle to send an electrical current to the area. This current causes rhythmic muscle contractions. They are two non-invasive physical therapy technologies that may help stimulate muscle contractions or treat pain. Both are sometimes called pulse stimulation.

## EMS has the following major benefits



## TENS is used to relieve pain and muscle spasms caused by a wide range of conditions including:



# RLT

Red & infrared light therapy (RLT) is also a non-invasive technology which uses specific wavelengths in the red (normally 660nm) & infrared light(850nm) spectrum to treat a variety of conditions in muscle, bone & skin. Red light can be absorbed into the skin to a depth of about eight to 10 millimeters, at which point it has positive effects on cellular energy and multiple nervous system and metabolic processes. It's been scientifically approved to have the following benefits.



We combine RLT technology with EMS and TENS to develop some innovative belt stimulators which are mainly for double effective pain relief, stretch marks removal, skin rejuvenation and toning muscles. More new products with RLT technology will be released soon in coming months.



## Heated Ab Toning Belt Stimulator

Item No.: S21-B19B

Target Area: Waist and Abdomen

Applicable Waist Size:

Approximately, 23.7-43.7"/60-111cm

### Intended Use:



MUSCLE  
TONING



BOOST BLOOD  
CIRCULATION



HEATING



MASSAGE

### Product Features:

1. 2 in 1, EMS+Heating ensure better toning, strengthening and firming abdominal muscles with greater results.
2. Skin friendly and breathable PU with Lycra material brings you good wrap and strong support.
3. Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
4. Wireless control through the USB rechargeable console and RF remote controller.
5. 3-level temperature adjustment with precise NTC temperature control and real-time temperature monitoring to avoid low-temperature burns.
6. 12 different modes with 15 adjustable intensity levels available.
7. It can be used for 8 days continuously if one time(25 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse + Heating
Main material	PU+Lycra
Working mode	12 modes
Pulse intensity	15 levels
Pulse frequency	0-120Hz
Heating	3 levels
Heating temperature	39/41/43 °C ±2
Automatic shutdown	About 25 minutes
Usage time	About 200 minutes
Battery Type	Controller: Rechargeable polymer lithium battery Remote control: button battery
Battery capacity	2,000mAh
Charging port	Type-C
Product Size	Belt: about 950*120mm Extension strap: about 350*65mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× Remote control
- 1× USB Charing Cable
- 1× User Manual





## Ab Toning Belt Stimulator for postpartum moms

Item No.: S28B-B27

Target Area: Abdomen

Applicable Waist Size:

Approximately, 28.8-53.2"/73-135cm

### Intended Use:



FAT REDUCTION



MUSCLE REBUILDING



STRETCH MARKS REMOVAL



PAIN RELIEF



RED LIGHT THERAPY



BOOST BLOOD CIRCULATION

### Product Features:


- 2 in 1, Pluse+RLT features ensure abdominal muscle toning, stretch marks reduction, pain relief, wound healing and blood circulation increase.
- Skin friendly and breathable Lycra material brings you good wrap and strong support.
- Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
- Handheld controller with USB rechargeable battery inside and digital display, easy to control and use it discreetly.
- 42 pcs LED units (660nm RED lighth+ 850nm Infrared) with total 126 pcs wicks.
- 7 different modes with 15 adjustable intensity levels available.
- It can be used for 5 days continuously if one time(25 minutes auto-shutdown) in a day after being fully charged.
- Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse + Red Light Therapy
Main material	ABS+Lycra
Working mode	7 modes
Pulse intensity	15 levels
Pulse frequency	0-100Hz
Red Light Wavelength	660nm*2+850nm*1 42 units/126 wicks
Red light level	1 level
Red light temperature	43°C ± 2
Automatic shutdown	About 25 minutes
Usage time	About 110 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	2,500mAh
Charging port	Type-C
Product Size	Belt: about 1,100*173mm Extension strap: about 400*70mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× USB Charging Cable
- 1× User Manual





## EMS Heating Belt for Ab and Lower Back

Item No.: S16-B3

Target Area: Waist and Abdomen

Applicable Waist Size:

Approximately, 24-44"/60-114cm

### Intended Use:



MUSCLE  
TONING



BOOST BLOOD  
CIRCULATION



HEATING



MASSAGE

### Product Features:

1. 2 in 1, EMS+Heating ensure better toning, strengthening and firming abdominal muscles with greater results.
2. Skin friendly and breathable Lycra material brings you good wrap and strong support.
3. Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
4. Wireless control through the USB rechargeable battery inside and digital display, easy to control and use it discreetly.
5. 3-level temperature adjustment with precise NTC temperature control and real-time temperature monitoring to avoid low-temperature burns.
6. 7 different modes with 15 adjustable intensity levels available.
7. It can be used for 8 days continuously if one time(25 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse + Heating
Main material	ABS+Lycra
Working mode	7 modes
Pulse intensity	15 levels
Pulse frequency	0-100Hz
Heating	3 levels
Heating temperature	39/41/43 °C ±2
Automatic shutdown	About 25 minutes
Usage time	About 200 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	2,000mAh
Charging port	Type-C
Product Size	Belt: about 950*122mm Extension strap: about 350*75mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× USB Charging Cable
- 1× User Manual





## Full Arm EMS Exerciser

Item No.: A8

Target Area: Full Arms

Main Users: Adults

### Intended Use:



MUSCLE  
TONING



BOOST BLOOD  
CIRCULATION



BUILD MUSCLES  
WITHOUT  
HEAVY LOADS



STONE UP  
ANYTIME  
ANYWHERE

### Product Features:

1. Building arm strength and endurance, which helps you to achieve better arm exercise effect no need heavy loads.
2. Massage and relax tired muscles, massaging full arms after doing workouts.
3. Improve blood circulation to prevent joint injuries and no gel replacement needed.
4. Wireless control through the USB rechargeable battery inside and digital display, easy to control and use it discreetly.
5. Light-weighted but solid, small size, very portable.
6. 5 different modes with 15 adjustable intensity levels available.
7. It can be used for 14 days continuously if one time (15 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, after exercising, and when traveling.

### Specifications

Function	Pulse
Main material	ABS
Working mode	5 modes
Pulse intensity	15 levels
Pulse frequency	0- 140Hz
Automatic shutdown	About 15 minutes
Usage time	About 210 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	350mAh
Charging port	Micro USB
Product Size	About 168*135*39mm

### WHAT' S IN THE BOX?

- 1× Product Unit
- 1× USB Charing Cable
- 1× User Manual





## Wearable EMS Face Slimming Device

Item No.: F3

Target Area: Face

Main Users: Adults

### Intended Use:



TIGHTEN  
AND LIFT  
YOUR FACE



ACHIEVE A  
PERFECT FACIAL  
CONTOUR



BOOST BLOOD  
CIRCULATION



TONE UP  
ANYTIME  
ANYWHERE

### Product Features:

1. Tighten and lift your face to achieve a perfect facial contour and improve the elasticity of the upper and lower jaw skin.
2. 8 stainless floating electrodes offer great skin contact on your cheeks.
3. High quality elastic silicon rubber band ensures good stability on head and easy to put it on.
4. 3 different modes with 10 adjustable intensity levels available.
5. USB Type-C rechargeable battery.
6. It can be used for 11 days continuously if one time(10 minutes auto-shutdown) in a day after being fully charged.
7. More comfortable to use it with a wet facial mask.
8. Easy to wipe and clean by wet cloth.

### Specifications

Function	Pulse
Main material	ABS+Silicone+Stainless
Working mode	3 modes
Pulse intensity	10 levels
Pulse frequency	0-200Hz
Automatic shutdown	About 10 minutes
Usage time	About 110 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	180mAh
Charging port	Type-C
Product Size	Net unit (fully open): 180*180*47mm Headband : 555*51*3mm

### WHAT'S IN THE BOX?

- 1× Product Unit
- 1× USB Charging Cable
- 1× User Manual







## 4 in 1 Waterproof Body Contour Beauty Device

Item No.: F9

Target Area: Whole Body

Main Users: Woman

### Intended Use:



LIFT YOUR CHEST



MUSCLE REBUILDING



STRETCH MARKS REMOVAL



IMPROVE SKIN



RED LIGHT THERAPY



BOOST BLOOD CIRCULATION

### Product Features:


1. 4 in 1 function combining CA+RF+EMS+RLT technology to tone and firm body skin while smoothing away cellulite and wrinkles.
2. Can be used on whole body, especially for breast health through the BR mode.
3. IPX7 waterproof feature ensures you to enjoy bathtime without worry.
4. 3 different modes with 5 adjustable intensity levels available.
5. USB Type-C rechargeable battery.
6. It can be used for 18 days continuously if one time(10 minutes auto-shutdown) in a day after being fully charged.
7. Long handle design makes you to reach your back easily.
8. Can be rinsed under tap water.

### Specifications

Function	CA + RF + EMS + RLT
Main material	ABS+Stainless
Working mode	3 modes
Pulse intensity	5 levels
IPX	7 levels
CA frequency	3MHz
RF frequency	1MHz
Pulse frequency	0-50Hz
Automatic shutdown	About 10 minutes
Red light Wavelength	620nm/625nm-15 units
Usage time	30-180 minutes (depending on different modes)
Battery Type	Rechargeable polymer lithium battery
Battery capacity	1,200mAh
Charging port	Type-C
Product Size	about 217*78*49mm (unit only)

### WHAT'S IN THE BOX?

- 1 × Product Unit
- 1 × USB Charging Cable
- 1 × User Manual





## Wireless TENS Pain Relief Mini Pad Unit

Item No.: S3-P2A-2U

Target Area: Whole Body

Main Users: Adults

### Intended Use:



MUSCLE  
TONING



BOOST BLOOD  
CIRCULATION



MASSAGE



MUSCLE  
REBUILDING

### Product Features:


1. The FDA Cleared and ANVISA Certified TENS mini pad unit can be used anywhere on your body for muscle and joint pain relief including Arthritis, osteoarthritis, tendonitis etc and also helpful for increasing muscle training and blood circulation.
2. Medical grade and skin friendly gel ensures good contact without irritation and allergy.
3. 2 different modes with 15 adjustable intensity levels available.
4. Wireless controller with USB rechargeable battery.
5. It can be used for 7 days continuously if one time(25 minutes auto-shutdown) in a day after being fully charged.
6. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse
Main material	ABS+PU
Working mode	2 modes
Pulse intensity	15 levels
Pulse frequency	0-41Hz
Automatic shutdown	About 25 minutes
Usage time	About 180 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	90mAh
Charging port	Micro USB
Product Size	About 200*70mm

### WHAT'S IN THE BOX?

- 1× Control Unit
- 1× Pad With Gels
- 1× USB Charging Cable
- 1× User Manual





## 3 IN 1 Lower Back Pain Relief Brace Belt Stimulator

Item No.: S28A-B21D

Target Area: Lower Back  
 Applicable Waist Size:  
 Approximately, 29.6-56.7"/75-144cm

### Intended Use:



PAIN RELIEF



SUPPORT  
LUMBAR SPINE



HEATING



BOOST BLOOD  
CIRCULATION

### Product Features:

1. 3 in 1, Pluse+Heating+Tourmaline features relieve muscle stiffness and joint pain, increase muscle endurance, better support lumbar spine and boost blood circulation.
2. Skin friendly and breathable Lycra material brings you good wrap and strong support.
3. Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
4. Handheld controller with USB rechargeable battery inside and digital display, easy to control and use it discreetly.
5. 3-level temperature adjustment with precise NTC temperature control and real-time temperature monitoring to avoid low-temperature burns.
6. 7 different modes with 15 adjustable intensity levels available.
7. It can be used for 9 days continuously if one time (25 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse + Heating + Tourmaline
Main material	ABS+Lycra
Working mode	7 modes
Pulse intensity	15 levels
Pulse frequency	0-100Hz
Tourmaline	42 pcs (tourmaline)
Heating	3 levels
Heating temperature	39/41/43 °C ±2
Automatic shutdown	About 25 minutes
Usage time	About 220 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	2,500mAh
Charging port	Type-C
Product Size	Belt: about 1,100*215mm Extension strap: about 400*100mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× USB Charging Cable
- 1× User Manual





## Heated Lower Back Brace Belt Stimulator

Item No.: S28A-B16B

Target Area: Lower Back  
 Applicable Waist Size:  
 Approximately, 24.8-49.2"/63-125cm

### Intended Use:



PAIN RELIEF



SUPPORT  
LUMBAR SPINE



MASSAGE



BOOST BLOOD  
CIRCULATION

### Product Features:

1. Heating with pulse features can relieve muscle stiffness and joint pain, increase muscle endurance, better support lumbar spine and boost blood circulation.
2. Skin friendly and breathable Lycra material brings you good wrap and strong support.
3. Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
4. Handheld controller with USB rechargeable battery inside and digital display, easy to control and use it discreetly.
5. 3-level temperature adjustment with precise NTC temperature control and real-time temperature monitoring to avoid low-temperature burns.
6. 7 different modes with 15 adjustable intensity levels available.
7. It can be used for 10 days continuously if one time (25 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

### Specifications

Function	Pulse + Heating
Main material	ABS+Lycra
Working mode	7 modes
Pulse intensity	15 levels
Pulse frequency	0-100Hz
Heating	3 levels
Heating temperature	39/41/43 °C ±2
Automatic shutdown	About 25 minutes
Usage time	About 250 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	2,500mAh
Charging port	Type-C
Product Size	Belt: about 1,000*190mm Extension strap: about 400*70mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× USB Charging Cable
- 1× User Manual





## Red Light Therapy Brace Belt for Whole Body

Item No.: S30A-B22

Target Area: Whole Body  
 Main Users: Adults  
 Applicable Waist Size:  
 Approximately, 23.7-59.1"/60-150cm

### Intended Use:



PAIN RELIEF



IMPROVE SKIN



RED LIGHT THERAPY



REDUCE INFLAMMATION


### Product Features:

1. Red light therapy belt is great for muscle and joint pain relief, wrinkles removal, inflammation reduction, wound healing and blood circulation.
2. Skin friendly and breathable Lycra material brings you good wrap and strong support.
3. 108 pcs LED units (660nm RED light + 850nm Infrared light) with total 324 pcs wicks.
4. 2-way power supply optional, through a wire controller connecting a power bank /AC adaptor OR a controller with 10,000mAh rechargeable battery inside.
5. 2 different modes with 3 adjustable levels of temperature setting and 5 levels of timer setting(10/15/20/25/30mins).
6. It can be used for 16 days continuously if one time (10 minutes auto-shutdown) in a day after being fully charged.
7. Easy to wipe and clean with wet cloth.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling.

Specifications	
Function	Red light and near-infrared light
Main material	ABS+Lycra
Working mode	2 modes Mode 1: Always on Mode 2: 10Hz frequency flashing
Red Light Wavelength	660nm*2+850nm*1 108 units/324 wicks
Red light timer	5 levels 10/15/20/25/30 minutes
Red light temperature	3 levels 39/41/43 °C±2
Usage time	10 minutes by default after power on about 160 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	10,000mAh
Charging port	Type-C
Product Size	Belt: about 1,157*170mm Extension strap: about 450*170mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× USB Charging Cable
- 1× User Manual





## 2 IN 1 Lower Back Brace Belt Stimulator with Red Light Therapy and EMS/TENS function

Item No.: S30B-B26

Target Area: Lower Back

Applicable Waist Size:

Approximately, 30.7-58.3"/78-148cm

### Intended Use:



PAIN RELIEF



SUPPORT LUMBAR SPINE



IMPROVE SKIN



RED LIGHT THERAPY



REDUCE INFLAMMATION



BOOST BLOOD CIRCULATION

### Product Features:


1. 2 in 1, red light therapy with pulse features can relieve muscle stiffness and joint pain, increase muscle endurance, better support lumbar spine and boost blood circulation.
2. Skin friendly and breathable Lycra material brings you good wrap and strong support.
3. 105 pcs LED units (660nm RED light + 850nm Infrared light) with total 315 pcs wicks.
4. Soft and comfortable electrodes fully fit the skin with perfect contact and no gel replacement needed.
5. Handheld controller with USB rechargeable battery inside and digital display, easy to control and use it discreetly.
6. 7 different modes with 15 adjustable intensity levels available.
7. It can be used for 16 days continuously if one time (10 minutes auto-shutdown) in a day after being fully charged.
8. Feel comfortable to use it anytime, anywhere at home, in the office, during sports, and when traveling

### Specifications

Function	Pulse + Red Light Therapy
Main material	ABS+Lycra
Working mode	EMS: 7 modes TENS: 3 modes RLT: 1 modes
Pulse intensity	15 levels
Pulse frequency	0-100Hz
Red Light Wavelength	660nm*2+850nm*1 105 units/315 wicks
Timer	5 levels 10/15/20/25/30 minutes
Red light level	1 level
Red light temperature	43°C±2
Usage time	10 minutes by default after power on about 160 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	10,000mAh
Charging port	Type-C
Product Size	Belt: about 1,138*220mm Extension strap: about 450*100mm

### WHAT'S IN THE BOX?

- 1× Main Device
- 1× Belt
- 1× Extension Belt
- 1× USB Charging Cable
- 1× User Manual





## EMS Foot Massager

Item No.: S3R-T01

Target Area: Lower limbs

Main Users: Adults

### Intended Use:



PAIN RELIEF



MUSCLE TONING



FAT REDUCTION



BOOST BLOOD CIRCULATION

### Product Features:

1. The FDA Cleared and ANVISA Certified foot massager can be used feet and leg for muscle and joint pain relief including Arthritis, osteoarthritis, tendonitis etc and also helpful for increasing muscle training and blood circulation.
2. 2 different modes with 15 adjustable intensity levels available.
3. Wireless control through the USB rechargeable console and RF remote controller.
4. It can be used for 7 days continuously if one time(25 minutes auto-shutdown) in a day after being fully charged.
5. The foot mat is Non-adhesive, silver layer grids as conductive electrodes, very easy to be cleaned.
6. Feel comfortable to use it anytime, anywhere at home, in the office, after exercising, and when traveling.

Specifications	
Function	Pulse
Main material	ABS+PU
Working mode	2 modes
Pulse intensity	15 levels
Pulse frequency	0-41Hz
Automatic shutdown	About 25 minutes
Usage time	About 180 minutes
Battery Type	Controller: Rechargeable polymer lithium battery Remote control: button battery
Battery capacity	90mAh
Charging port	Micro USB
Product Size	About 310*310*17mm

### WHAT'S IN THE BOX?

- 1× Main Unit
- 1× Remote control
- 1× Foot Mat
- 1× USB Charing Cable
- 1× User Manual





## EMS Colorful Foot Mat

Item No.: S11

Target Area: Lower limbs

Main Users: Adults

### Intended Use:



PAIN RELIEF



FAT REDUCTION



MASSAGE



BOOST BLOOD CIRCULATION


### Product Features:

1. The colorful foot mat can be used feet and leg for muscle and joint pain relief including Arthritis, osteoarthritis, tendonitis etc and also helpful for increasing muscle training and blood circulation.
2. It can also training your arms when your hands put on it.
3. 8 different modes with 15 adjustable intensity levels available.
4. Wireless controller with USB rechargeable battery.
5. Digital display and touch screen, easy to control and use it discreetly.
6. It can be used for 12 days continuously if one time(15 minutes auto-shutdown) in a day after being fully charged.
7. The foot mat is Non-adhesive, silver layer as conductive electrodes, soft & Rollable & Super slim, very easy to be cleaned.
8. Feel comfortable to use it anytime, anywhere at home, in the office, after exercising, and when traveling.

Specifications	
Function	Pulse
Main material	ABS+TPU
Working mode	8 modes
Pulse intensity	15 levels
Pulse frequency	0-70Hz
Automatic shutdown	About 15 minutes
Usage time	About 180 minutes
Battery Type	Rechargeable polymer lithium battery
Battery capacity	300mAh
Charging port	Type-C
Product Size	about 360*338*19mm

### WHAT' S IN THE BOX?

- 1× Main Device
- 1× USB Charing Cable
- 1× Fabric Pouch
- 1× User Manual







## Maternity Belt

Item No.: B29

Target Area: Waist and Abdomen

Main Users: Pregnant Women

Applicable Waist Size:

S: 70-100cm/27.6"-39.4"

M: 80-100cm/31.5"-39.4"

L: 90-120cm/35.4"-47.2"

XL: 100-130cm/39.4"-51.2"

XXL: 110-140cm/43.3"-55.1"

### Intended Use:

Maternity belts are worn by pregnant women to make their pregnancies more comfortable. These pregnancy supports can also reduce pain and discomfort due to conditions like hip separation, hip instability, lordosis and lymphedema.

### Product Features:

1. Wraps around your belly and waist, lessening the burden of pregnancy weight.
2. Helps you get relief from hip, pelvic and sacroiliac pain, while further minimizing spinal strain.
3. High elastic fabric is breathable, soft, comfortable and safe, and can be adjusted at will.
4. The portable hook and loop fastener is strong and durable.
5. Multiple colors and sizes available to choose.
6. Suitable for use in the middle and late stages of pregnancy.



#### Specifications

Material	Polyester + elastic fiber
Color	Beige/black/pink/white



## 3 in 1 Postpartum Support Belt

Item No.: B30

Target Area: Waist and Abdomen

Main Users: Pregnant Women

Applicable Waist Size:

M: 74-84cm/29.1"-33.1"

L: 85-94cm/33.5"- 37"

XL: 95-104cm/37.4"-40.9"

XXL: 105-114cm/41.3-44.9"

### Intended Use:

Postpartum belts are supportive garments that are worn around the waist after giving birth. These are designed to help reduce swelling, provide support to the lower back and abdominal muscles, and promote healing after delivery.

### Product Features:

1. 3-in-1 wrapping support function, including waist belt, belly belt, and pelvic belt.
2. Relieve back pain, correct posture, and tighten the pelvis.
3. High elastic fabric is breathable, soft, comfortable and safe, and can be adjusted at will.
4. The portable hook and loop fastener is strong and durable.
5. Multiple colors and sizes available.
6. Suitable for leisure and entertainment, housekeeping, shopping and other scenarios.

#### Specifications

Material	Polyester + elastic fiber
Color	Beige/black





## Back Support Belt with Removable Waist Pad

Item No.: B31

Target Area: Lower Back

Applicable Waist Size:

S: 71-80cm/28"-31.5"

M: 81-90cm/31.9"-35.4"

L: 91-100cm/35.8"-39.4"

XL: 101-115cm/39.8"-45.3"

XXL: 111-120cm/43.7"-47.2"

### Intended Use:

Back support belts can take pressure off the spine and lower back and make it easier to move without pain by providing compression, comfortable support, and preventing unnatural motions.


### Product Features:

1. Give strong support to your waist with removable waist pad.
2. Highly elastic and breathable fabric offers double pressure support.
3. The portable hook and loop fastener is strong and durable.
4. Lightweight and soft, can be worn under clothes.
5. Multiple colors and sizes available.
6. It can be used while doing yoga, weightlifting, running, housekeeping, and gardening.



#### Specifications

Material	Nylon + elastic band + PP support strip
Color	Black



## Back Support Belt with Metal Stays

Item No.: B32

Target Area: Lower Back

Applicable Waist Size:

S: 60-68cm/23.6"-26.8"

M: 70-82cm/27.6"-32.3"

L: 84-92cm/33.1"-36.2"

XL: 94-102cm/37"-40.2"

XXL: 104-110cm/40.9"-43.3"

### Intended Use:

Back Support Belts with Metal Stays adopt metal bionic support bars, combined with lightweight air mesh fabric, firmly support the spine for weightlifting, heavy lifting, relieve arthritis, disc herniation, sciatica, spinal stenosis, scoliosis, ankylosis, spondylitis, Back pain from surgery, spondylolysis, osteoporosis, etc.

### Product Features:

1. Give strong support to your waist with those built-in flexible metal stays.
2. Air mesh fabric and anatomical support strips will provide the back with high elasticity, breathability and durability to fit the curve of the lower back, ensuring a snug fit for the lower back, dynamic body balance and coordination of complex trunk movements.
3. The portable hook and loop fastener is strong and durable.
4. Lightweight and soft, can be worn under clothes.
5. Multiple colors and sizes available.
6. Suitable for people who sit and stand for a long time, for leisure and entertainment, and for doing housework.

#### Specifications

Material	Mesh cloth + hook and loop cloth + metal stays
Color	Gray/Black

