



H-4199.5F

Soil Grinder



Description:

Prepares soil samples to designated particle size for accurate, repeatable test results. Grinds one-pint sample in 15 seconds. Stainless steel construction. Includes a No. 10 perforated stainless plate. Shipping wt. 32 lbs. (15kg)

Supports the following standards:
ASTM D5102, ASTM D558, ASTM D698, ASTM D1557

Technical Specifications:





Model type: Lab Bench Oven, Digital
Max. Temp: 450°F (232°C),
Capacity: 10.6 cu. ft. (300L)
Power: 1920 watts, 230V 50/60Hz
Dimensions: Inside: 25.5" x 24" x 30" (648 x 610 x 762mm)
Overall: 33" x 35.5" x 34" (838 x 902 x 864mm)

Available models:

H-4199 Soil Grinder, 120V 60Hz
H-4199.5F Soil Grinder, 220V 50Hz



Accessories

Part Number & Name	
	H-4199.A No. 4 Perforated, Stainless Plate
	H-4199.B No. 10 Perforated, Stainless Plate
	H-4199.C No. 35 Perforated, Stainless Plate
	H-4199.D 2mm Perforated, Stainless Plate

Replacement Parts

Part Number & Name	
	H-4199.7 Beater Assembly