

H-4325.5F

Humboldt, Economy Sieve Shaker



Description:

The Humboldt economy sieve shaker can be used with 3", 5" and 8" sieves. It can handle up to ten 8" sieves, twelve 5" sieves, sixteen 3" full-height sieves or eighteen half-height 8" sieves. This shaker uses a 1/4 hp motor with a 30-minute timer. Unit should be bolted to bench for correct operation.

Technical Specifications:




Model: H-4325.5F
Power: 230V 50Hz
Dimensions: 15"w x 15"d x 45"h (380 x 380 x 1143mm).

Manufacturer:

Humboldt Mfg. Co.
Chicago, IL, USA



Product Options

Part Number & Name	
	H-4325 Humboldt, Economy Sieve Shaker, Economy 120V 60Hz
	H-4325.2F Humboldt, Economy Sieve Shaker, Economy 220V 60Hz
	H-4325.5F Humboldt, Economy Sieve Shaker, Economy 220V 50Hz