



# Humboldt Sieve Shakers

Manufactured by Humboldt Mfg. Co. - Chicago, IL, USA



## Description:

The Humboldt Sieve Shakers are built upon the economical and time-proven design, which has been providing dependable service for many years.

The economy sieve shakers feature a 30-minute, easy-to-use, mechanical timer and durable construction.

The digital version allows the user to set a time of up to 60 minutes with automatic shutoff with an accuracy of 0.5%.

They can be used with 3", 5", 8", 10" and 12" sieves. See the different models below for details.

Unit should be bolted to bench for correct operation.

## Model / Technical Specifications:

Suitable models for Australia:

- H-4325.5F Economy Sieve Shaker**  
 Operated via: 30 min mechanical timer  
 Capacity: 10 x 8", 12 x 5" or 16 x 3" full-height sieves  
 Dimensions: 15"w x 15"d x 45"h (380 x 380 x 1143mm).  
 Power: 230V 50Hz – ¼ hp motor
- H-4330.5F Large, Economy Sieve Shaker**  
 Operated via: 30 min mechanical timer  
 Capacity: 11 x 8", 7 x 10" or 7 x 12" full-height sieves  
 Dimensions: 21"w x 18"d x 47"h (533 x 4570 x 1194mm).  
 Power: 230V 50Hz – ¼ hp motor
- HA-4325V.5F Variable-Speed Sieve Shaker, Digital Timer**  
 Operated via: 60min digital timer  
 Capacity: 10 x 8", 12 x 5" or 16 x 3" full-height sieves  
 Dimensions: 15"w x 15"d x 45"h (380 x 380 x 1143mm).  
 Power: 230V 50Hz – 1/3 hp motor

- HA-4325.5F Sieve Shaker, Digital Timer**  
 Operated via: 60 min digital timer  
 Capacity: 10x8", 12x5" or 16x3" full-height sieves  
 Dimensions: 15"w x 15"d x 45"h (380x380x1143mm).  
 Power: 230V 50Hz – ¼ hp motor
- HA-4330.5F Large Sieve Shaker, Digital Timer**  
 Operated via: 60min digital timer  
 Capacity: 11x8", 7x10" or 7x12" full-height sieves  
 Dimensions: 21"w x 18"d x 47"h (533x4570x1194mm).  
 Power: 230V 50Hz – ¼ hp motor
- HA-4330V.5F Large Variable-Speed Shaker, D.Timer**  
 Operated via: 60min digital timer  
 Capacity: 11x8", 7x10" or 7x12" full-height sieves  
 Dimensions: 15"w x 15"d x 45"h (380x380x1143mm).  
 Power: 230V 50Hz – 1/3 hp motor



H-4325.5F



HA-4325.5F



HA-4325V.5F



H-4330.5F



## Hand-Operated Model / Technical Specifications:

### Humboldt, Hand-Operated Sieve Shaker ASTM C136

The Humboldt, hand-operated sieve shaker is ideal for locations where electricity is not available or limited. can be used with 3", 5" and 8" sieves. It can handle up to ten 8" sieves, twelve 5" sieves, sixteen 3" full-height sieves or eighteen half-height 8" sieves.


Unit should be bolted to bench for correct operation.

<b>H-4310</b>	<b>Humboldt, Hand-Operated Sieve Shaker</b>
Operated via:	Handle
Capacity:	10 x 8", 12 x 5" or 16 x 3" full-height sieves
Dimensions:	15"w x 15"d x 45"h (380 x 380 x 1143mm).
Power:	N/A



## Alternative Models:

There are several additional brands and models suitable for Australia. See our website for more details.



H-4324 Sound Enclosure, Ro-Tap



H-4326.5F Sieve Shaker for 3" Sieves



H-4321.5F Dura Tap™ Sieve Shaker



H-4328.5F Wet/Dry Sieve Shaker