





## PLUG PUSHERS



### TWO FLANGE PUSHER

The Adept Two Flange Plug Pusher is designed for fast and easy insertion of the Adept Two Flange and Soluble Two Flange plug.

- > Durable, lightweight and ergonomic
- > Polished finish for ease of cleaning and sterilization
- > Depth stop feature for accurate insertion and prevents operator contact with animal

Product Code	0142-060
Supplied	Single unit
Material	Stainless Steel / Plastic
Colour	White (additional colours available on request)
Dimensions	210mm
Plugs	Two Flange (1324-001) Soluble Two Flange (1324-002)

### PLUG 42 PUSHER

The Adept Plug 42 Pusher is designed for fast and easy insertion of the Adept Plug 42 plug.

- > Durable, lightweight and ergonomic
- > Polished finish for ease of cleaning and sterilization

Product Code	0142-020
Supplied	Single Unit
Material	Stainless Steel / Plastic
Colour	White (additional colours available on request)
Dimensions	210mm
Plugs	Plug 42 (0072-002)