


EVERYTHING YOU NEED TO KNOW!



PRICES



Large Tip
(Set of 2)
\$80.00



Small Tip
(Set of 2)
\$80.00



AS200E
\$15,999.00

SPECIFICATIONS

Handpiece	Weight	245g (with battery)	Tip	Size	22 x 18 mm (L) 18 x 16 mm (S)
	Dimensions	248 x 48 x 37 mm		Autoclavable	Up to 60 times 134°C 4mins
Scanning	Full-arch Accuracy	20 µm	Battery	Scan Duration	120 min
	Dimensions	16 x 14 mm (L) 12 x 12 mm (S)		Charging	Wireless
	Capture Button	Duo buttons	Calibration	Calibration free	



Our partnership with Voxel and NEOLab in the US and Canada is a testament to our commitment of pushing the boundaries of innovation to introduce new technologies. Together, we bring a wireless intraoral scanner optimized for pediatric dentistry with the goal of providing streamlined digital solutions across various clinical needs.



Voxel Dental is a dedicated team of dental technology consultants and support specialists, offering digital imaging and 3D printing solutions to dental practices across the United States. Our products enhance patient care while streamlining dental practice workflow. Digitally transform your practice with Voxel Dental Solutions.



First 3 Appliances for Free

With a long-standing history in pediatric appliances, NEOLab Junior is the technical laboratory partner you've been looking for. We're technology and quality driven, enabling us to bring you the most cutting-edge workflows with unmatched fit. Save time in your practice with our innovative solutions.



The *NEW* AS200E Wireless Intraoral Scanner



Optimized for Pediatric Dentistry
Calibration Free. Pristine Accuracy. Rapid Turnaround.

DIGITAL SCANNING TECHNOLOGY OPTIMIZED FOR PEDIATRIC DENTISTRY

Small in size, **HUGE IN POTENTIAL**

**5 YEAR
WARRANTY**

**NO
MONTHLY
FEES!**

**FREE
ON-BOARDING
& SUPPORT**

**CALIBRATION
FREE**

**DROP
PROTECTION**

Goop Free with Wireless Digital Scans

- Full arch scanning in 20 seconds
- Skip the alginates with a wireless high resolution intraoral scanner
- AI-powered enhancement
- True-to-life color

Expanding Platform

- Digitalize your workflow with labs you already use
- Access your data anywhere using our HIPPA compliant cloud services
- Online 3D viewer

Breakthrough Wireless Technology

- First wireless charging scanner on the market
- Up to 2-hours continuous scanning (equivalent to 60 cases in a single charge)
- Tip in two sizes with a small tip option optimized for pediatric use
- Sleek and lightweight design to ensure minimal hand and wrist strain

Seamless Integration

- Instant and easy transfer of scans with integrated dental lab (NEOLab Junior or your lab of choice)
- Easily track your orders through Scan Pro
- Track historical data on patients for ease of continuity and diagnosis

FIRST 3 APPLIANCES FOR FREE

When you choose **neolab**
junior



Watch Video

Band & Loop




Watch Video

Lingual Arch



Watch Video

Expander



Watch Video

KD Overlay

3D Printed Metal Appliances

NEOLab is a leader in metal printed appliances within patent-pending digital design process

- Designed to avoid both subgingival and interproximal interference
- Skip clinical appointments and avoid spacers
- 3D printed bands provide a better fit and stronger hold compared to traditional bands

Rx Wizards

The easiest way to submit pediatric prescriptions

- **Batch shipping** – group shipments together to save on shipping costs
- **Integrated tracking** – fully transparent tracking from each case from time of scan to arrival
- **Custom templates** – save your favorite designs to be used again and again
- **On demand shipping labels** – Sending a physical impression? Rx Wizard can generate as many shipping labels as you need with a single click