



9th Nov 2022

Dear customers,

HSA grants 24 months shelf life for Flowflex ART

On the 9th of November 2022, the Health Sciences Authority (HSA) has granted Flowflex an extension of its shelf life of 24 months from the manufacturing date.


This letter is to notify you that the Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing) kits have a longer shelf life than the expiry date stated on the packaging

The extension of 24 months expiry date applies to the following batches.

COV2015140; COV2015218; COV2015251; COV2015252; COV2015253; COV2015254; COV2015255; COV2015261; COV2015279; COV2015280; COV2025052; COV2025053; COV2025054; COV2025078; COV2025079; COV2025080; COV2025081; COV2025082; COV2025083; COV2025085; COV2025086; COV2025087; COV2025088; COV2025098; COV2025099; COV2025100; COV2025101; COV2025123; COV2025125; COV1110232; COV1120153.

How to determine the new expiry date:

Step 1: Refer to the manufacturing date printed on the cassette pouch.



Step 2: Add 24 months to the manufacturing date.

If the manufacturing date is 2022-02-15, the new expiry date is 2024-02-14

For more information, please visit our website at www.flowflex.sg