

## NEW!! INFECTION CONTROL



500ml or 50ml

### One For All Foaming Hand Sanitizer

Alcohol sanitizers don't protect you. One For All works differently: It cleanses the skin by killing 99.9% of germs in seconds and forms a protective barrier that keeps on killing germs until you wash it off.

### AiroPur Aerosol Suction Cup

Attaches to the hose of your HVE unit to control aerosols in the dental operator. It's simple to use, inexpensive enough to be disposable... and designed so it is fully autoclavable for multiple uses. Available in packs of 40 units.



40-ct

## VACU-MIXER



### Transfer Jig

Allows hands-free transfer of impression material from one cartridge to another.



2-ct

### Vacu-Mixer Auto Mixing Machine

Hand-held, portable and re-chargeable. Weighing less than 2 lbs, the Vacu-Mixer holds a 120-ml cartridge (rather than 50 ml). The larger cartridge lowers the cost of materials, leaves less waste, and delivers a perfect mix every time.

### Dynamic Mixing Tips

Saves material usually wasted in larger mixing tips. Less waste means bigger impression material cost savings.



50-ct

## VACCU-SIL IMPRESSION MATERIAL

In an extensive study performed by a Utah-based independent dental products evaluation organization, Vaccu-sil was rated number one in the world amongst all VPS impression materials for value and performance.



### Vaccu-sil 50-ml Cartridges

Heavy Regular or Fast | Light Regular or Fast

4 x 50 ml

## MOJO SYRINGE



### The Mojo Syringe (ready to fill)

The impression material wasted in a typical mixing tip is 2.4ml. That's enough to do 3-4 more crowns with our specially-designed Mojo Syringe. Plus, the extra-narrow tip allows you to apply material at the sulcus—circumferentially—without having to lift the tip, minimizing air bubbles and costly re-takes. Ergonomically designed for easy access to the distals of second molars. 100-ct.



### The Mojo Complete (pre-filled)

Our Mojo Syringe, prefilled with Vaccu-sil light-body VPS impression material. Eliminate cross-contamination and messy clean-ups. Light Body Fast or Regular 50-ct.



### Extra Mojo Tips

Handy when enough material remains in a Mojo Syringe for another impression. Eliminates wasted material. 10-ct.

### 120ml Cartridges for the Vacu-Mixer

Heavy Body Regular

Heavy Body Fast

10 x 120 ml



## HARMONY DUAL ARCH TRAYS



50-ct

### Harmony Posterior Dual Arch Tray

The tray's narrow, 2.6 mm stainless steel bar prevents tuberosity impingement and gives the tray rigidity equal to or better than any other disposable posterior dual arch trays.



40-ct

### Harmony Anterior Dual Arch Tray

It's an anterior arch tray + facial plane relator. Captures the prepared teeth, the opposing arch, the bite registration, and the patient's facial midline in one simple step. Can be used without alignment sticks for single units, or with them for multiple units. No more slanted crowns or veneers!

## HOETCH ACID ETCH GEL

A 37% phosphoric acid etch gel for use on dentin and enamel. It provides an excellent prepared dental surface for accepting resin at a fraction of the cost.



3 x 5ml | 30 x 5ml

## VACCU-SIL BITE REGISTRATION MATERIAL

2 x 50 ml



Like our top-rated Vaccu-sil VPS impression material, our Vaccu-sil Bite Registration stands up to the major brands at the lowest price.

## MYLAR MATRIX ROLL

300 feet



At just nine microns wide, our Mylar is ultra-thin. It can be used during cementation of ceramic restorations and removed afterwards without leaving a gap.

## VACUCORD GINGIVAL RETRACTION CORD 100 inches



A knitted, non-medicated retraction cord used for rapid tissue displacement and absorption of fluids at the sulcus, in preparation for the dental impression. Sizes 00, 0, 1 and 2, as well as a kit containing one each of all four.

## THE FACIAL PLANE RELATOR



Easily records for your laboratory technician the accurate bite registration and the exact position of the desired midline in relation to a horizontal plane in just 60 seconds. 12-ct | 24-ct

## SAND TRAP

Our patented clear-plastic device is designed to confine and suction away aluminum oxide particulates during intra-oral or chairside sandblasting. 12-ct | 100-ct



## 4 Reasons to Choose Ho Dental Company:

- ✓ Products created by a practicing dentist for dentists
- ✓ Up to 30% lower cost than competing brands
- ✓ Recognized by key industry evaluators for value and performance
- ✓ 35 years of product innovation