ENDO-MATE TC2 Cordless Endodontic Micromotor

Smarter & Safer

Work Smarter and Safer with Freedom of Movement

A large LCD, simple 5-key operation and a lightweight, cordless handpiece assure easy use even during the most delicate endodontic procedures. The Auto Reverse & Alarm Function alerts you with an audible sound to let you know that the load is about to reach the preset torque level allowing you to unload the file even before the Auto Reverse is setting in. 5 programs for different file systems. TC2 supports most major brands of Ni-Ti files.

Mini head & Slim neck Miniature head allows access even to the most difficult to reach areas. Slim neck assures higher visibility. Stress free for more comfort Cordless handpiece delivers freedom of movement. Lightweight, ergonomic design. LCD panel Large LCD panel and access-friendly key.

Program function

User definable 5 presets - Speed, Torque, Gear Ratio & Auto Reverse mode. Last program settings can be memorized when the power turns on.







Safety for delicate endodontic treatment

Smart Auto Reverse & Alarm Function for safer endodontic procedures. The feedback circuit keeps stable rotation even when the load on the motor handpiece changes. Power turns off automatically to improve safety and save energy. Calibration function for more precise rotation speed and torque setting.

This function cannot completely absorb the individual difference of the actual

Suitable for use with most major brands of Ni-Ti files

Selecting gear ratio attachment applies the entire rotation and torque.

Rotation Speed (min ⁻¹)	Gear Ratio	Torque Level (Ncm)	
200-900	10 : 1	0.2-2.2	
140-550	16:1	0.3-3.0	
100-450	20:1	0.4-4.0	

Two kinds of alarm sound for safer and efficient treatment

Alarm enables operators to prevent to put too much burden on files, realizing less metal fatigue.

When load reaches the preset torque, the rotating file stops and Auto reverse mode activates.			MAX	Auto Reverse On !
At about 75% of the preset torqu alarm activates.	ie,	75%		Alarm
At about 50% of the preset torque, alarm activates.	50%			Alarm
) Preset	t Torque Leve		100	

ENDO-MATE TC2 Complete Set

ENDO-MATE TC2 Complete MPA Set

MODEL	ORDER CODE
ENDO-MATE TC2 MPA	Y1001028

MPA-F16R Head • ENDO-MATE TC2 Motor • Battery Charger

- Specifications
 Speed range: 2,000-9,000 min⁻¹
- Charge: 1.5 h (depending on workload)

 Continuous use: 1.2 h (depending on workload)

