

# Lace Up Ankle Brace with Figure 8 Strapping

SKU: AKB177

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



## Instructions For Use

**Step 1.** Insert foot and tighten laces.

**Step 2.** Wrap ankle with velcro straps.

**Step 3.** Secure velcro straps for secure fit.



Still need help?  
**Perfect Fit Guarantee!**

Visit our website or scan the QR code to find out more.

### Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse thoroughly and Air Dry.
- Do NOT dry clean, iron, or bleach.  

