

Low Back Support Belt

SKU: BKB122

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Open your back brace completely. Locate the Velcro closure on the waist panel.

Step 2. Attach the end of the extension belt to the velcro closure.

Step 3. Wrap around your waist and secure the waist panels with Velcro.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Spray with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.   

