

Patella Knee Strap

SKU: KSP197

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Unfasten Velcro and open the brace.

Step 2. Wrap the strap around your knee, just below the knee cap. The padded portion should sit where the patellar tendon attaches to the patella.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Spray with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

