

Pneumatic Foot Device

SKU: PFD210

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Place foot into strap.

Step 2. Wrap strap and secure with velcro.

Step 3. If needed, trim strap to desired length.

*Ensure proper fit prior to trimming. Once trimmed, product can not be returned.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

