

Mini Thumb Spica Splint

SKU: WRW258

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Unfasten all straps place thumb into orthosis.

Step 2. Velcro the thumb strap so it fits securely.



Step 3. Velcro the wrist strap so it fits securely.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  



BraceDirect.com • 6860 Big Ridge Rd • Hixson, TN 37343 • (800) 257-5870