

Breezy Basix Wheelchair

Lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look

- ✓ Available in both self propelled and transit models
- ✓ Different size options to choose from
- ✓ Fixed and half-folding backrest available
- ✓ Armrests conveniently flip back, for easier transfers
- ✓ Legrests - swing away, height adjustable, removable
- ✓ Footplates – plastic flip up with heel loops, angle adjustable
- ✓ Swing away and adjustable Anti-Tips
- ✓ Rear wheels, solid, quick release
- ✓ Australian standards and crash tested approved



Breezy Basix Transit Wheelchair with Half Folding Back Rest

Product Code	074104-001	074104-003	074104-005	074104-006
Colour	Silver	Silver	Silver	Silver
Overall Height (cm)	72	72	72	72
Overall Depth (cm)	94	94	94	94
Overall Width (cm)	58	63	68	73
Seat Width (cm)	38	43	48	52
Seat Depth (cm)	41-49	41-49	41-49	41-49
Seat Height (cm)	43	43	43	43
Armrest Height (cm)	20.5-23	20.5-23	20.5-23	20.5-23
Backrest Height (cm)	37-42	37-42	37-42	37-42
Max Load (kg)	125	125	125	125
Net Weight (kg)	14.9	14.9	14.9	14.9
Warranty (months)	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]

Breezy Basix Transit Wheelchair with Fixed Back Rest

Product Code	074129-001	074129-003	074129-005	074129-006
Colour	Silver	Silver	Silver	Silver
Overall Height (cm)	72	72	72	72
Overall Depth (cm)	94	94	94	94
Overall Width (cm)	58	63	68	73
Seat Width (cm)	38	43	48	52
Seat Depth (cm)	41-46	41-46	41-46	41-46
Seat Height (cm)	46-51	46-51	46-51	46-51
Armrest Height (cm)	20.5-23	20.5-23	20.5-23	20.5-23
Backrest Height (cm)	41-46	41-46	41-46	41-46
Max Load (kg)	125	125	125	125
Net Weight (kg)	14.9	14.9	14.9	14.9
Warranty (months)	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]

Breezy Basix Manual Wheelchair with Half Folding Back Rest




Product Code	074103-001	074103-003	074103-005	074103-006
Colour	Silver	Silver	Silver	Silver
Overall Height (cm)	72	72	72	72
Overall Depth (cm)	94	94	94	94
Overall Width (cm)	58	63	68	73
Seat Width (cm)	38	43	48	52
Seat Depth (cm)	41-46	41-46	41-46	41-46
Seat Height (cm)	46-51	46-51	46-51	46-51
Armrest Height (cm)	20.5-23	20.5-23	20.5-23	20.5-23
Backrest Height (cm)	41-46	41-46	41-46	41-46
Max Load (kg)	125	125	125	125
Net Weight (kg)	14.9	14.9	14.9	14.9
Warranty (months)	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]


Breezy Basix Manual Wheelchair with Fixed Back Rest

Product Code	074101-001	074101-003	074101-005	074101-006
Colour	Silver	Silver	Silver	Silver
Overall Height (cm)	72	72	72	72
Overall Depth (cm)	94	94	94	94
Overall Width (cm)	58	63	68	73
Seat Width (cm)	38	43	48	52
Seat Depth (cm)	41-46	41-46	41-46	41-46
Seat Height (cm)	46-51	46-51	46-51	46-51
Armrest Height (cm)	20.5-23	20.5-23	20.5-23	20.5-23
Backrest Height (cm)	41-46	41-46	41-46	41-46
Max Load (kg)	125	125	125	125
Net Weight (kg)	14.9	14.9	14.9	14.9
Warranty (months)	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]	60 [12 on castors, fixings, seat, straps, and wheels]

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au