

## Configura Advance Chair

Highly configurable specialist chair. Innovative design for complex and challenging seating needs.



- ✓ Tool-less adjustments
- ✓ Removable, angle and height-adjustable footplate
- ✓ 37 degrees of tilt in space and full length compensating legrest
- ✓ 130 degrees backrest recline, 10 degrees of forward tilt
- ✓ Visco integrated self-supply cushion
- ✓ Lockable, reduced shear backrest recline to maximise pelvic stability

### Extras:

- ✓ Upgradeable backrest options
- ✓ Drop-down armrest option (helps with transfers and sling insertion)
- ✓ Forward tilt option (standard on electric chair)
- ✓ **Manual or electric option**



Product Code	CR5435	CR5438
Mechanism	<b>Manual Tilt-in-Space</b>	<b>Electric Tilt-in-Space</b>
Colour	Charcoal	
Fabric	Vinyl fabric with VP on all areas	
Seat Width Options	350mm (14in) – 550mm (22in)	
Seat Depth Options	400mm (16in) – 550mm (22in)	
Seat Height Options	330mm (13in) – 530mm (21in)	
Weight Capacity (SWL)	160kg	

We care about  
your comfort!

[www.capeabilities.com.au](http://www.capeabilities.com.au)

📍 26 Strickland Street, Bunbury WA 6230  
📞 08 9717 6622  
✉ [bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

📍 Unit 1/64 West Street, Busselton WA 6280  
📞 08 9751 1446  
✉ [sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)