

## Configura Comfort Bariatric Lift Recliner Chair

Brings Configura technology to bariatric seating (Twin motor recliner chair designed to provide comfortable seating for users with mild to moderate postural requirements)



### Three sizes and 3 colours to choose from

- ✓ Increase in measurements for comfortable seating
- ✓ Duratek fabric provides improved pressure reduction
- ✓ Pressure reducing fabric on backrest and armrest as standard
- ✓ Size adjustment options available including depth and height
- ✓ Interchangeable backrest options – adjustable pillow backrests
- ✓ Pressure care seating options available – including ROHO cushions
- ✓ 2-button handset & carer control options available
- ✓ Standard & Tilt In Space options

Product Code	Standard			Tilt In Space		
Sizes	16"	18" 254KG	18" 317KG	16"	18"	20"
Colour/Code	<b>CR5412</b>	<b>CR5414</b>	<b>CR5420</b>	<b>CR5410</b>	<b>CR5416</b>	<b>CR5418</b>
Seat Height (cm)	40.6	45.7	45.7	40.6	45.7	50.8
Seat Width (cm)	60.9/66/71.1	60.9/66/71.1	76.2/ 81.2/ 86.3	60.9/66/71.1	60.9/66/71.1	60.9/66/71.1
Seat Depth (cm)	45.7/50.8/ 55.8	45.7/50.8/ 55.8	45.7/50.8/ 55.9	45.7/50.8/ 55.8	45.7/50.8/ 55.8	45.7/50.8/ 55.8
Armrest Height (cm)	20	20	20	20	20	20
Overall Height (cm)	115	115	115	115	115	115
Overall Width (cm)	105	105	105	105	105	105
<b>Max Load (kg)</b>	<b>254</b>	<b>254</b>	<b>317</b>	<b>254</b>	<b>254</b>	<b>254</b>
<b>Leg-rest weight capacity (kg)</b>	<b>127</b>	<b>127</b>	<b>158</b>	<b>127</b>	<b>127</b>	<b>127</b>
Warranty (years)	2	2	2	2	2	2

We care about  
your comfort!

[www.capeabilities.com.au](http://www.capeabilities.com.au)

26 Strickland Street, Bunbury WA 6230  
08 9717 6622  
[bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

Unit 1/64 West Street, Busselton WA 6280  
08 9751 1446  
[sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)