

Configura Comfort Lift Recliner Chair

Twin motor recliner chair designed to provide comfortable seating for users with mild to moderate postural requirements

Three sizes and 3 colours to choose from




- ✓ Pressure reducing fabric on backrest and armrest as standard
- ✓ Size adjustment options available including depth and height
- ✓ Interchangeable backrest options – adjustable pillow backrests
- ✓ Pressure care seating options available – including ROHO cushions
- ✓ 2-button handset & carer control options available
- ✓ Tilt in space as standard



Product Code	Small	Medium	Large
Colours/Fabric	Vinyl-Black (CR5443) Duratek (CR5404) Taupe (CR5422)	Vinyl-Black (CR5444) Duratek (CR5406) Taupe (CR5424)	Vinyl-Black (CR5445) Taupe (CR5426)
Seat Height (cm)	38.8 / 40.6 / 43.1 / 45.7	43.9 / 45.7 / 48.2 / 50.8	43.9 / 45.7 / 48.2 / 50.8
Seat Width (cm)	45.7	50.8	55.9
Seat Depth (cm)	40.6 / 45.7 / 50.8	40.6 / 45.7 / 50.8	40.6 / 45.7 / 50.8
Armrest Height (cm)	18	18	18
Overall Height (cm)	109.2	114.3	114.3
Overall Width (cm)	71.1	76.2	81.2
Overall Length - Upright Position (cm)	94	94	94
Overall Length - Recline Position (cm)	170.2	170.2	170.2
Chair to Wall/ Furniture (cm)	45.7	48.3	48.3
Max Load (kg)	130	160	160
Warranty (years)	2	2	2

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au