

Cura1 Falls Prevention Products – Cordless Range

Cura1 is a leading manufacturer of these products, offering 12-month warranty and in-house customer and technical support. Electronic communications has informed the innovative range of cura1 technology.




Monitor (Controller/Receiver) + Cordless Pad Transmitter + Cordless Item

The cordless range is a simple operation and easy to clean. The 50m range of the monitor allows monitoring of up to 6 devices (including chair pad, bed pad, wheelchair chair pad, floor mat, and Activecare pager). These items are used for carers within the same home, supplying innovative technology to prevent falls.

Other Cura1 product ranges include: (Speak to a member of staff to explore options)

- 2558 Wireless Call Button with Reset & Pull Cord
- 3011 SeizureAlert Monitor with Alarm Relay Output (inc. pager option)
- 2900 DoorAlert + Door Monitor



We care about
your comfort!

www.capeabilities.com.au

26 Strickland Street, Bunbury WA 6230
08 9717 6622
bunbury@capeabilities.com.au

Unit 1/64 West Street, Busselton WA 6280
08 9751 1446
sales@capeabilities.com.au

Item	Product Code	Measurements	Information
Six Channel Cordless Falls Monitor	2700	120 x 77 x 40mm	<ul style="list-style-type: none"> Monitors up to 6 wireless devices. Includes batteries, durable plastic enclosure & caregiver reset key. Amber LED visual alarm indicator. Up to 6 months battery life. 3 selectable alarm tones & adjustable alarm volume. "Pause" feature to allow caregiver to attend to the patient without setting the alarm off. Can be connected to a nurse system.
Three Channel Cordless Falls Monitor	2702	120 x 83 x 35mm	<ul style="list-style-type: none"> Monitors up to 3 cordless devices. Amber LED visual alarm indicator. Up to 6 months battery life. 3 selectable alarm tones & adjustable alarm volume. Can be connected to a nurse system.
Cordless Pad Transmitter	2719	433MHz operating frequency & 50m range	<ul style="list-style-type: none"> Simple operation and maintenance (12-month battery life). Fits inside cordless pads.
Cordless Chair Pad	2630	380 x 255mm	<ul style="list-style-type: none"> Maximum weight capacity: 200kg Easy to clean, made from non-allergenic soft PVC and has vent holes in the neck to reduce the risk of moisture ingress. For use with Cura1 2719 Transmitter.
Cordless Bed Pad	2629	760 x 510mm	<ul style="list-style-type: none"> Maximum weight capacity: 200kg Easy to clean, made from non-allergenic soft PVC and has vent holes in the neck to reduce the risk of moisture ingress. For use with Cura1 2719 Transmitter.
Cordless Under Mattress Pressure Pad	3515	800 x 290mm	<ul style="list-style-type: none"> Maximum weight capacity: 200kg Under mattress means completely non-intrusive bed monitoring. Easy to clean, made from non-allergenic soft PVC and has vent holes in the neck to reduce the risk of moisture ingress. For use with Cura1 2719 Transmitter.
Cordless Fall Mat	2624	1220 x 610mm	<ul style="list-style-type: none"> Maximum weight capacity: 200kg Easy to clean, made from heavy duty non-slip material (surface material is foam PU (Polyurethane) and bottom material is PVC without Phthalate). For use with Cura1 2719 Transmitter.